

Mimic pathogen

Aim of the experiment

The purpose of this experiment is to express viral nucleic acid in *Escherichia coli* to simulate the state that viral nucleic acid is wrapped in a protein shell under real conditions. We constructed a plasmid with a virus sequence (about 100bp). It was introduced into *E. coli* and induced with IPTG to produce target RNA for subsequent experiments.

(Note: The virus plasmid only contains a small segment of the virus sequence and is not pathogenic. It is safe and legal.)

Materials

- *E. coli* BL21 competent cells
- LB medium
- 500ml flasks
- Kanachloramphenicol (Km, 100ug/ml) : Dissolve 0.1g of Km by adding it to 8 mL of deionized water and then vortexing. Add deionized water to bring the volume to 10 mL and filter-sterilize with a 0.22-μm syringe filter. Store at -20 °C.
- Shaker

Procedure

1. Thaw one vial of *E. coli* BL21 competent cells on ice for 30 min, and then add 1 μL of 50 ng/μL of H1N1 plasmid. Incubate on ice for 5 min.
 2. Heat-shock the cells by placing the vial into a 42 °C pre-heated water bath for 45 s, and then cold-shock the cells on ice for 2 min.
 3. Add 200 μL of LB medium to the cells and plate 100 μL of cell suspension on a pre-warmed
-

LB plate containing 100 µg/mL Kanachloramphenicol. Incubate the plate overnight in a 37 °C incubator.

4. The next day, inoculate 5 mL of LB medium containing two antibiotics with a single colony and incubate the culture overnight at 37 °C in a biological shaker at 300 r.p.m.
 5. Inoculate 500 mL of LB medium, containing two antibiotics and determine the optical density (OD, 600 nm). Shake cultures at 37 °C, 300 r.p.m.
 6. Monitor the OD every hour until the cells reach an OD of 0.4–0.6, and then induce expression by adding 1 mL/L 0.5 M IPTG and shake the cultures overnight at 300 r.p.m. in apre-chilled 28°C biological shaker to make E. coli transcribe RNA.
-