



UPDATED 20.10.2020 18:47



by Franka Butzbach (created by Philip Schulz)



Introduction

Digestion Ligation is a procedure in restriction cloning and it is the procedure with which we assemble the constructs in our iGEM project. ALWAYS: when we do DIGLIG, which needs 4h time then first check if one of the small blocks A or B is available! In this protocol, we use the insertion of MFa into the p02 Vector as an example

Materials

- › MQ
- › 10X Tango Buffer
- › ATP
- › T4 DNA Ligase
- › Enzymes
- › BSA
- › T4 Buffer (as an alternative to 10X Tango Buffer and ATP)
- › Plasmids
 - › Backbone and insert

Procedure

How to Calculate the Concentrations in the case of our MFalpha:

1. <https://nebiocalculator.neb.com/#!/ligation>

Go to this website and type in the DNA concentrations and bp lengths that you have for insert and vector, **choose 2:1 ratio**. So in this example we need around 13ng of insert.

Note: For LII assembly, You have to put in the bp length of the whole plasmid and not that of the insert it contains.

Ligation Calculator

This tool will calculate the mass of insert required at several molar insert:vector ratios in the range needed for typical ligation reactions.

Insert DNA length	Required insert DNA mass
<input style="width: 80%;" type="text" value="317"/> bp	6.401 ng (1:1)
<input style="width: 80%;" type="text" value="4952"/> bp	12.80 ng (2:1)
<input style="width: 80%;" type="text" value="100"/> ng	19.20 ng (3:1)
	32.01 ng (5:1)
	44.81 ng (7:1)

2. The length of our Mfa is 317bp (including the amplification overhangs) and we choose (2:1) ratio, so you need 12.8ng of insert Mfa.

3. Nanodrop your insert: we have 68.5ng/μL Mfa

Calculate how much of the stock DNA solution needs to be added. in the case of Mfa:

$(12.8\text{ng}) / (68.5\text{ng}/\mu\text{L}) = 0.19\mu\text{L}$. This is too small a volume to pipette directly. If this happens, you need to dilute it (for example, 1.9μL of the insert in 8.1μL of MQ) This would yield a solution with Mfa concentration of 12.8ng/μL so you could just add 1μL of this dilution to your reaction to add 12.8ng.

Use the formula $c_1 \times V_1 = c_2 \times V_2$ to dilute your inserts.

V_1 = amount of stock you add to dilution

C_1 = conc. of stock (NANODROPED)

V_2 = final volume = added stock plus MQ to dilute (10μL)

C_2 = conc. of dilution (=required insert DNA mass)

$V_1 \mu\text{L} + (10\mu\text{L} - V_1 \mu\text{L}) \mu\text{LMQ}$ = dilution (total volume of 10μL)

LI construct DIGLIG

	A	B	Example for MFa in p02
1	H2O	fill up to 15 ul	10.5 µl
2	Vector (7.5 nM)	calculate as described above	dilute 1 ul in 3 ul MQ, then add 1 ul
3	Insert (15 nM)	calculate as described above	dilute 1.9 ul in 8.1 ul MQ, then add 1 ul
4	10X Tango Buffer	1 µl	1 µl
5	T4 DNA Ligase (5 units/µl)	0.25 µl	0.25 µl
6	Enzyme (5 units/µl)	0.25 µl	0.25 µl
7	ATP (10 mM)	1 µl	1 µl

LII constructs DIGLIG

	A	B
1	H2O	fill up to 15 ul
2	Vector (7.5 nM) p(10) for all!	calculate as described above
3	Insert A-B Promoter (15 nM)	calculate as described above
4	Insert B-C MFa or dummy BB6 for B+ (15 nM)	calculate as described above
5	Insert C-D Receptor (15 nM)	calculate as described above
6	Insert D-E Fluorescent protein(15 nM) 125.9 ng (2:1)	calculate as described above
7	Insert E-F Terminator (15 nM)	calculate as described above
8	Insert F-G BB9 dummy (15 nM)	calculate as described above
9	10X Tango Buffer	1 .5 µl

10	T4 DNA Ligase (5 units/ μ l) our tube has 400000 u per ml \rightarrow 400u/uL	0.25uL
11	Enzyme (5 units/ μ l)	0.25 μ l dilution: (1uL EZ in 3uL MQ \rightarrow to volume 4uL) \rightarrow then add 1uL
12	ATP (10 mM)	1 μ l

4. The program on the thermocycler is **diglig long**, which can also be used for short insert as well as for the LII constructs.

