

Experiment: Improve's Transformation: Banana-amilCP, pp21-amilCP, original-amilCP
 (~~Banana~~CP) (pp21CP) (OriCP)
 Ban CP

Resuspend's Concentration:

	260/230	260/280	Conc. ng/ul
original-amilCP	4.16	1.62	33.29
Banana-amilCP	3.46	1.62	24.31
pp21-amilCP	1.78	1.24	70.42

Dilute to 5 ng/ul

original-amilCP: 2 ul + 10 ul ddH₂O to 5ng/ul

Banana-amilCP: 2 ul + 8 ul " " "

pp21-amilCP: 1 ul + 13 ul " " "

Transformation:

1. Add ³³~~55~~ ul of E.coli-DH5-α competent cells into an eppendorf.
2. Add 3.5 ul of targeted plasmid into the eppendorf.
3. Cool the eppendorf on ice for 3-5 min.
4. Heatshock the eppendorf at 42°C for 1 min.
5. Cool the eppendorf on ice for 3-5 min.
6. Add 500 ul of antibiotics-free LB media into the eppendorf.
7. Incubate the eppendorf at 37°C for 1 hour. (With shaking)
8. Centrifuge the eppendorf at 15000 rpm for 2 min.
9. Remove 436 ul of the supernatant of the eppendorf.
10. Vortex and spin down.

11. Take 100 μ l of the transformed bacteria and add it on a LB plate with Kanamycin
12. Incubate the plate for overnight. (8:32 pm $\rightarrow$ 10:32pm)

* For step number 11, we made the following modifications with suggestions from Phil

- We add competent cells into 6 dishes:

- ① pp21_AmilCP 5 μ l (less than 5 μ l)

- ② pp21_AmilCP 95 μ l

- ③ Banana_AmilCP 5 μ l

- ④ Banana_AmilCP 95 μ l

- ⑤ Original_AmilCP 5 μ l

- ⑥ Original_AmilCP 95 μ l

- Reason: Since this is a retransformation of IDT Product, there ~~is~~ are possibilities that our transformed competent cells will be too prolific after overnight incubation if we add 100 μ l of competent cells on to a petri dish. So we add 2 different amount of competent cells, 5 μ l and 95 μ l, to insure that the bacteria has a lower chance to outgrowth the petri-dish

SUBJECT

DATE

2020.08.21 楊永毅

Place: Ming Dao High school

Temp: 25°C

Pressure: 1 atm

Experiment: Culture Bacteria of

- ori CP $\odot$ & $\ominus$
- ban CP $\odot$ & $\ominus$
- pp21 CP $\odot$ & $\ominus$

Protocol:

~~1. Add 7 ml LB + Kan into 6 centrifuge tubes~~

~~2. Add~~

1. follow 2020.05.25 culturing bacteria protocol.

SUBJECT

DATE

2020.08.22 Shen

Place: Ming Dao High School

Experiment: ① Bacteria storage ② Plasmid extraction

① 300ul of 50% glycerol
+ 600ul of O/N culture bacteria

→ -80°C Third floor → IGEN 2020 CSNU BOX

Item: ori CP ori CP ban CP ban CP PP21 CP PP21 CP
x6 0822 0822 0822 0822 0822 0822
① ② ① ② ① ②

②	260/280	260/230	Conc. ng/ul		260/280	260/230	Conc. ng/ul
ori CP ①	1.47 2.89	3.89	195.67	ori CP ②	1.60	4.63	141.52
ban CP ①	1.57	4.44	151.27	ban CP ②	1.56	4.34	150.57
PP21 ①	1.40	3.79	214.38	PP21 ②	1.46	4.10	191.49

Place: Ming Dao High School

Experiment: Improve of AmilCP (Banana Control [BanCP0], Original [OriCP0])

Plate B ~~State~~

	J1	J2	J3
	Banana Control	Banana Exp	Original
DNase Free	0.97 μ L 2 μL	0.6 μ L	1.04 μ L
Sol A	2 μ L	2 μ L	2 μ L
Sol B	1.5 μ L	1.5 μ L	1.5 μ L
RNase Inh.	0.2 μ L	0.2 μ L	0.2 μ L
Plasmid	0.33 μ L	0.33 μ L	0.26 μ L
Trigger	0	0.37 μ L	0

* Trigger conc. 132.57 ng/ μ L, 260/230: 4.76, 260/280: 1.68

Plate B

	J5	J6	J7	J8	J9
	Ban_Control	Banana_Exp	Original	RFP_Control	Water
DNase Free	1.94 μ L	1.19 μ L	2.09 μ L	0.94 μ L	10 μ L
Sol A	4 μ L	4 μ L	4 μ L	2 μ L *	0
Sol B	3.15 μ L *	3 μ L *	3 μ L *	1.5 μ L	0
RNase Inh.	0.4 μ L *	0.4 μ L *	0.4 μ L *	0.2 μ L *	0
Plasmid	0.66 μ L	0.66 μ L	0.51 μ L	0.36 μ L	0
Trigger	0 μ L	0.75 μ L	0 μ L	0 μ L	0

* RFP₀ conc.: 140.82 ng/ μ L, 260/230: 3.76, 260/280: 1.94

Results: Failed, Please refer to 2020.08.24 Improve absorbance on the drive.

Discussions:

- ① We exclude the interfere of synthesis kit by measuring the excitation and emission spectrum of mRFP, which turned out to be normal.
- ② We suspected that the issues may stem out from the process of maturation for amilCP after protein synthesis.

2020.08.25 Ly

Place: Ming Dao High School

Temp: 25°C

Pressure: 1 atm

Experiment: Improve of Am-CP

- ① original (ori-CP)
- ② improved (ban-CP) + trigger
- ③ improved (ban-CP)

Plate B	ori-CP		ban-CP+tri		ban-CP		H ₂ O		
	J ₁₁	J ₁₂	J ₁₃	J ₁₄	J ₁₅	J ₁₆	J ₁₇	J ₁₈	J ₁₉
DNase free H ₂ O	1.04	4.6	0.6	2.4	0.97	3.88	5	0	0
A	2	8	2	8	2	8	0	0	0
B	1.5	6	1.5	6	1.5	6	0	0	0
RNase inhibitor	0.2*	0.8	0.2*	0.8	0.2*	0.8	0	0	0
Plasmid	0.26	1.04	0.33	1.32	0.33	1.32	0	0	0
Trigger	0	0	0.37	1.48	0	0	0	0	0
total	5	20	5	20	5	20	5	20	0

SUBJECT

DATE

楊承睿·潘建安
2020.08.26 Ivy·張桓睿

Place: Ming Das High School

Temp: 25°C

Pressure: 1 atm

Experiment: Improve of Ami/CP: ① ban CP + 1x tri
② ban CP + 2x tri
③ ban CP without tri

Protocol:

1. Plate B

	J21	J22	J23	J24
DNase free H ₂ O	2.4	0.92	3.88	20
A	8	8	8	0
B	6	6	6	0
RNase inhibitor	0.8	0.8	0.8	0
Plasmid	1.32	1.32	1.32	0
Trigger	1.48	2.96	0	0
total	20	20	20	20

2. Plate reader

wavelength: 588nm

Temp: 37°C

Time: 4:00:00, Interval: 00:02:00

Orbital: 0:01:00

Result: Increase trigger ~~CAN'T~~ CAN'T increase the yield of CP.

SUBJECT

DATE

2020.08.26 楊承睿·張程睿

Place: Ming Dao High School

Temp: 25°C

Pressure: 1atm

Experiment: ~~Improve of~~ Culture bacteria of Banana Trigger 1 & 2

Protocol:

1. Add 7 ml of LB media with ampicillin into a centrifuge tube.
2. Take out frozen bacteria from -80°C and place it on ice.
3. Use a tooth pick to pick the bacteria, and mix thoroughly with the LB media.
4. Incubate overnight.

SUBJECT

DATE

張桓睿·陳中霖

2020.08.27 陳虹好

Place: Mingdao High School

Temperature: 24°C

Experiment: ~~For~~ Plasmid Extraction of Banana Trigger; Triple Replication of Improve

0827 BT1

0827 BT2

* Banana Trigger 1: 166.98 ng/ul * Banana Trigger 2: 178.22 ng/ul

Plate B	K1	K2	K3	K4	K5
	Ban_Tri	Ban_Con	Ori	Water	Blank
DNase-Free	2.68 ✓	3.88 ✓	4.18 ✓	20 ✓	
Sol A	8 ✓	8 ✓	8 ✓	0	
Sol B	6 ✓	6 ✓	6 ✓	0	
RNase-Inh	0.8 ✓	0.8 ✓	0.8 ✓	0	
Plasmid	1.32 ✓	1.32 ✓	1.02 ✓	0	
Trigger	1.20 ✓	0	0	0	
Total	20ul	20ul	20ul	20	X
	K6	K7	K8	K9	Blank K10
DNase Free	2.68 ✓	3.88 ✓	4.18 ✓	20 ✓	
Sol A	8 ✓	8 ✓	8 ✓	0	
Sol B	6 ✓	6 ✓	6 ✓	0	
RNase-Inh	0.8 ✓	0.8 ✓	0.8 ✓	0	
Plasmid	1.32 ✓	1.32 ✓	1.02 ✓	0	
Trigger	1.20 ✓	0	0	0	X
	K11	K12	K13	K14	K15
DNase Free	2.68 ✓	3.88 ✓	4.18 ✓	20 ✓	
Sol A	8 ✓	8 ✓	8 ✓	0	
Sol B	6 ✓	6 ✓	6 ✓	0	
RNase-Inh	0.8 ✓	0.8 ✓	0.8 ✓	0	
Plasmid	1.32 ✓	1.32 ✓	1.02 ✓	0	