

Saliva Collection and RNA amplification

Reagent	Source	Identifier
QuickExtract DNA Extraction Solution	Lucigen	QE09050
NASBA liquid kit	Life Sciences Advanced Technologies Inc.	SKU: NWK-1
Tris (1 M), pH = 8, RNase free	Invitrogen	AM9855G
Sodium Hydroxide	Sigma-Aldrich	71687
1M MgCl ₂	Invitrogen	AM9530G
2M KCl	Invitrogen	AM9640G
DTT	Sigma-Aldrich	43816
DMSO	Sigma-Aldrich	276855
dNTP set, 100 mM	Invitrogen	10297018
NTP set, 100 mM	Thermo Scientific	R0481
Bio-11-UTP (75 mM)	Invitrogen	AM8451
RNase H	NEB	M0297L
ProtoScript II reverse transcriptase	NEB	M0368S
T7 RNA polymerase	NEB	M0251L
BSA 20 mg/ml	NEB	B9000S
Direct-zol RNA Miniprep Plus	Zymo Research	R2070
PCRD lateral flow assay	Abingdon Health	FG-FD51673
Qubit RNA HS Assay Kit	Invitrogen	Q32852
PowerUp™ SYBR™ Green qPCR Master Mix	Applied Biosystems	15340939

1. Lysis of saliva samples Mix crude saliva at 1:1 ratio with QuickExtract DNA Extraction Solution. Incubate at 95 °C for 5 min to ensure complete lysis of virus and inactivation of proteinase K.

2. NASBA reaction Take 1 µl from the product of Step 1 and add into the NASBA reaction mixture to make a total volume of 20 µl. Reaction mixture can either be prepared in-house or from the Life Sciences NASBA liquid kit (see tables below). a. Reaction mixture without the enzyme mix is incubated at 65 °C for 2 min followed by a 10-min incubation at 41 °C. Following that, 5 µl enzyme mix is added into the reaction and incubated at 41 °C for a further of 90-120 min. b. Alternatively, reaction mixture including the enzyme mix is prepared and incubated directly at 41 °C for a total of 90-120 min.

3. At the end of the NASBA reaction, RNA is purified from the end-product using Direct-zol RNA Miniprep kit, and eluted with 30 μ l of RNase free water. After purification, 4.2 μ l of purified RNA is mixed with 1.8 μ l of FAM labelled RNA capture oligo and 84 μ l of PCRD extraction buffer (see table below). Take 75 μ l of mix to the sample well of a PCRD test cassette. Results will be shown within 10 min.

Life Sciences reaction mixture (RM)

	vol.	stock conc.	conc. in RM
sample	1 μ l		
primers	1 μ l	0.5 μ M each	25 nM each
water +/- beacon	3 μ l		20 nM for beacon
buffer (NECB-24)	6.7 μ l		
nucleotide (NECN-24)	3.3 μ l		
enzyme mix (NEC-1-24)	5 μ l		
total volume	20 μ l		

In-house reaction mixture (RM)

	vol.	stock conc.	conc. in RM
sample	1 μ l		
primers	1 μ l	0.5 μ M each	25 nM each
water +/- beacon	4 μ l		20 nM for beacon
buffer with DMSO*	5 μ l		
nucleotide mix*	4 μ l		
enzyme mix*	5 μ l		
total volume	20 μ l		

References

1. Wu, Q., Suo, C., Brown, T., Wang, T., Teichmann, S. A., & Bassett, A. R. (2020). INSIGHT: a population scale COVID-19 testing strategy combining point-of-care diagnosis with centralised high-throughput sequencing. Cold Spring Harbor Laboratory.
<https://doi.org/10.1101/2020.06.01.127019>
2. Qin, Z., Peng, R., Baravik, I. K., & Liu, X. (2020). Fighting COVID-19: Integrated Micro- and Nanosystems for Viral Infection Diagnostics. *Matter*, 3(3), 628–651.
<https://doi.org/10.1016/j.matt.2020.06.015>