

PCR clean-up

Aim of the Experiment

This protocol can be used for the quick purification of an exponentially PCR amplified DNA of interest.

Materials

- PCR amplified DNA of interest
- Axygen PCR-clean up kit
- Eluent

Procedure

1. Adding 120 ul Buffer PCR-A into the product and centrifuged for 1 min with $12000 \times g$. Discarding the filtrate.
 2. 700ul Buffer W2 was added into the preparation tube and centrifuged for 1 min with $12000 \times g$. Discarding the filtrate.
 3. Repeat the above steps once and let stand for 5min (volatile ethanol).
 4. Add 30ul Eluent into the preparation tube, let stand for one minute, and centrifuge with $12000 \times g$ for 1min.
-