

Heat shock Transformation for TOP10 E.coli

Materials

- Chemically competent cells
- Purified plasmid
- Petri dish- with relevant selective antibiotic
- 5X KCM
 - 0.5M KCl
 - 0.15M CaCl₂
 - 0.25M MgCl₂

Procedure

1. Thaw Top10 *E.coli* cells on ice for 15 minutes
2. Add 2 µL KCM in a new 1.5 mL tube
3. Add 15-30 ng of plasmid (or 1-8 µL DNA mix or the ligation sample)
4. Add distilled water or ligation mix to a total volume of 10 µL
5. Incubate on ice for 2-5 minutes
6. Add 10 µL of the thawed competent cells.
7. Incubate on ice for 20 minutes
8. Heat shock at 42°C for 1 minute

9. Put on ice for 3 minutes
10. Add 180 μL of prewarmed TSB/LB or any nutritionally rich medium
11. Phenotype with shaking for 1 hour at 37°C
12. Plate on petri dish (160 μL on one half, 20 μL on the other half)
13. Incubate overnight at 37°C.