



UPDATED 9.10.2020 16:00



by Franka Butzbach (created by Philip Schulz)



Introduction

Use the Miniprep kit to extract DNA from liquid overnight cultures.

Materials

› MINIPREP KIT

- › and the RNase stored in 4°C fridge

› Additional materials

- › 1.5mL tubes, all pipettes, MQ, gloves

Procedure

Miniprep (omega biotek kit - quick guide)

<https://www.omegabiotek.com/wp-content/uploads/2017/08/D6945-Quick-Guide-011217.pdf>

1. grow 1-5mL culture overnight in a 10-20mL culture tube (should be happy shaking in the shaker incubator)
2. Centrifuge at 2,000 x g for 10 minutes at room temperature 22°C (Big centrifuge). Wear gloves:
Decant the liquid in the liquid Biowaste (needs to get autoclaved) -> we need white sandy bacteria stuff at the bottom
3. Add 250 µL Solution I mixed with RNase A (stored in fridge at 4°C). Vortex to mix thoroughly.
4. Transfer suspension into a new 1.5mL microcentrifuge tube.
5. Add 250µL Solution II. *Gently* invert the tube by hand several times to obtain a clear lysate, (gently so don't disrupt the DNA, avoid vigorous mixing. Do this for 2-3 min. DO NOT exceed 5 min!
6. Add 350µL Solution III. Immediately invert several times until a flocculent white precipitate forms.
Centrifuge at max speed (13,000xg) for 10 min. A compact white pellet will form.
7. Insert a HiBind[®] DNA Mini Column into a 2mL Collection Tube. Pay attention to not touch anything with your hands where there will be DNA afterwards, so don't contaminate stuff.
8. Transfer the cleared supernatant liquid from Step 5 by CAREFULLY pipetting it into the HiBind[®] DNA Mini Column+ 2mL tube (should be around 850µL). Pay attention to not touch or suck in the white pellet (this is just protein waste -> discard the tube in Biowaste autoclaved) Centrifuge at max speed (13'000) for 1 minute. Discard the liquid filtrate in the liquid Biowaste and reuse the collection tube.
The plasmid is in the column filter membrane.
9. Add 500 µL HBC Buffer (yellow bottle) diluted with 100% isopropanol Centrifuge at max speed for 1 min. Discard the filtrate and reuse the collection tube. NOTE: if there is a cross on the bottle cap it is already mixed with iso and ready to use
10. Add 700 µL DNA Wash Buffer (green bottle) diluted with 100% ethanol. Centrifuge at max speed for 1 min. Discard the filtrate and reuse the collection tube. NOTE: if there is a cross on the bottle cap it is already mixed with et and ready to use)
11. Centrifuge the empty HiBind[®] DNA Mini Column II at max speed for 2 min to dry the column. This step is critical for removal of trace ethanol that may interfere with downstream applications.
12. Transfer the HiBind[®] DNA Mini Column II into a nuclease-free 1.5mL microcentrifuge tube. (our normal tubes)
13. Add 30-100 (kyle: 50) µL Elution Buffer or sterile deionized water. Let sit at room temperature for 60 seconds. Centrifuge at max speed for 1 min. ATTENTION: the centrifuge goes anticlockwise so the caps of the tubes need to be in the correct direction otherwise they get destroyed by the force.

14. Store eluted DNA at -20°C in fridge