

Type IIS cloning (MoClo)

Introduction

Protocol adapted from "A Modular Cloning System for Standardized Assembly of Multigene Constructs" by Weber *et al.* by Dorottya Marko. Contributions from Björn Ancker Persson and Núria Garriga Alonso.

Materials

- 2uL Reaction mixture
 - T4 ligase buffer 10x 1 uL
 - BSA (1mg/mL) 0.5 uL
 - Bsal or SapI (5U/uL) 0.25 uL
 - T4 ligase (5U/uL) 0.25 uL
- DNA mixture
 - 15 fmol/uL insert (each insert) -> $\text{Length(bp)}/100 = \text{total ng needed for 10uL reaction}$
 - 7.5 fmol/uL -> $\text{Length(bp)}/200 = \text{total ng needed for 10uL reaction}$
- H₂O until 10uL

Procedure

Temperature Cycle

1. 37°C 3min Restriction
2. 16°C 4min Ligation
3. Go to step 1 x26
4. 50°C 5min
5. 80°C 10min Enzyme deactivation
6. 4°C Forever