

# Transformation of *E. coli*

## Introduction

Transformation of *E. coli* with USER Cloning assembled product (using competent cells).

## Materials

- Incubator at 37 °C
- Warm agar plates (media and supplements)
- USER reaction products
- 1.5 mL Eppendorf tubes

## Procedure

1. Aliquot 50 µL of chemically competent cells into 1.5 mL Eppendorf tubes (one for every transformation plus 2 for the positive and negative controls) on ice.
2. Warm selection plates to 37°C.
3. Add 10 µL of the chilled assembly product to the competent cells (or all the mixture). Mix gently by pipetting up and down or by flicking the tube 4-5 times. Do not vortex.
4. Place the mixture on ice for 30 minutes. Do not mix.
5. Heat shock at 42 °C for 30 seconds. Do not mix.
6. Transfer tubes to ice for 2 minutes.
7. Spread all the mixture (around 60 µL) of cells onto the selection plates.
8. Remember to make positive and negative controls.
9. Incubate overnight at 37 °C.