# Riskbedömning

## SB/IB Sterile work/LAF-bench

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#### Riskkategori

d: Låg risk;

The valve for the natural gas should be closed when not in use. The user must read through the SOP for the lamina flow hood with burner before using.

#### Exponeringsnivå

Låg;

#### Skyddsventilation

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Agens klass 1;

Working in LAF bench can protect sample from contamination and also avoid spores from spore forming fungi to be spread in the laboratory. If you are working with filamentous fungi, only use the dedicated LAF bench.

AFS 2005:1

Risk group 1 contains biological agents which do not normally cause infections in humans, together with non-pathogenic strains of pathogenic biological agents. Risk group 1 also includes biological agents which do not cause infection but can cause other ill-health, such as hypersensitivity or toxicogenicity, which is not connected with infections.

#### Farokoder

Natural gas  
Risk code R12 Highly flammable.  
The gas is supplied from a central system in the building.

### Lokaler

---

Ventilation; Annat;

### Personlig skyddsutrustning

---

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### Beskriv teknisk utrustning

---

SOP for the LAF bench and the burner is available next to the bench. Make sure to clean the bench with 70% Ethanol before and after your work. Remove your things from the bench (samples, pipette tips, sterile loops, etc.) once you are done.

### Avfallshantering

---

Radioaktivt avfall;

Liquid biological waste is collected in assigned container. Solid biological wastes are collected in assigned bags to be autoclaved and disposed as general waste.

Solid chemical wastes are discarded in the assigned box and liquid chemical wastes are collected in a container, then sent to a company for proper disposal.

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Natural gas

#### ACCIDENTAL RELEASE MEASURES

Use protective clothing. Wear self-contained breathing apparatus when entering area unless atmosphere is proved to be safe. Evacuate area. Ensure adequate air ventilation. Eliminate ignition sources. Try to stop release. Prevent from entering sewers, basements and workpits, or any place where its accumulation can be dangerous. Ventilate area.

Respiratory protection is not needed unless gas displaces the air and causes deficiency of oxygen and possibility of asphyxiation. Symptoms may include loss of mobility/consciousness. In low concentrations may cause narcotic effects. Symptoms may include dizziness, headache, nausea and loss of co-ordination. Remove victim to uncontaminated area wearing self contained breathing apparatus. Keep victim warm and rested. Call a doctor. Apply artificial respiration if breathing stopped. Incomplete combustion may form carbon monoxide

Do not extinguish a leaking gas flame unless absolutely necessary. Spontaneous/explosive re-ignition may occur. Extinguish any other fire.

Fire-fighting

Dry chemicals. Carbon dioxide (CO<sub>2</sub>). Water spray. Alcohol resistant foam.  
Improper extinguishing media Do not use water jet.

Inhalation; Fresh air and rest. Get medical attention if any discomfort continues. Perform artificial respiration if breathing has stopped.

Skin contact; Remove contaminated clothes and rinse skin thoroughly with water.

Eye contact; Rinse the eye with water immediately. Make sure to remove any contact lenses from the eyes before rinsing. Contact physician if discomfort continues.

Ingestion; Drink a few glasses of water or milk. Do not induce vomiting. Contact physician if larger quantity has been consumed.

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---

Ventilationsavbrott; Annat;

Keep gas/flame at proper distance to flammable chemicals.

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---

Use protective clothing. Localize (have in close vicinity) fire-extinguisher

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---

1. Acceptabel risk;

**Signatur Ansvarig chef**

Christer Larsson

**Datum**

---

**Riskbedömning : SB/IB Sterile work/LAF-bench**

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**2022-09-01**