

# Ni-NTA agarose protein purification

## Aim of the experiment

The purpose of this experiment is to purify the Cas13 protein, which has a 6×His tag and can be purified with Ni-NTA agarose. The whole experiment process should be carried out at 4°C to prevent protein denaturation.

## Materials

### instrument:

- Vortex oscillator
- Ultrasonic disruptor
- Rotating incubator QB-328
- Microplate reader

### Reagent:

- Phosphate Buffered Saline (1×) (PBS Buffer)
  - Phenylmethyl sulfonyl fluoride (PMSF)
  - Tris-HCl
  - NaCl
  - Tris(2-carboxyethyl)phosphine (TCEP)
  - Imidazole
  - Ultra-pure water
  - 20% ethanol
  - rTEV protease
  - SUMO protease
  - 50 mL column
  - Ni-NTA Agarose
  - Centrifugal filter
-

- Ultrafiltration tube
- BCA Protein Quantitation Assay Kit
- Buffer:

**Table 1: Washing Buffer A**

<b>Reagent</b>	<b>concentration</b>
Tris-HCl pH 8.0	50 mM
NaCl	30 mM
TCEP	1mM
Imidazole	10 mM
Constant volume to 1L, and degas ultrasonically	

**Table 2: Washing Buffer B**

<b>Reagent</b>	<b>concentration</b>
Tris-HCl pH 8.0	50 mM
NaCl	30 mM
TCEP	1mM
Imidazole	50 mM
Constant volume to 1L, and degas ultrasonically	

**Table 3: Elution Buffer**

<b>Reagent</b>	<b>concentration</b>
Tris-HCl pH 8.0	50 mM
NaCl	30 mM
TCEP	1mM
Imidazole	250 mM

---

Constant volume to 1L, and degas ultrasonically

---

**Table 4: Washing Buffer C**

Reagent	concentration
Tris-HCl pH 8.0	50 mM
NaCl	30 mM
TCEP	1mM
Imidazole	500 mM
Constant volume to 1L, and degas ultrasonically	

---

**Table 5: PB Buffer**

Reagent	concentration
Na <sub>2</sub> HPO <sub>4</sub> ·12H <sub>2</sub> O	50 mM
NaH <sub>2</sub> PO <sub>4</sub> ·2H <sub>2</sub> O	50 mM
Constant volume to 1L	

---

## Procedure

1. Mix the thawed protein thawed cell pellet with PBS buffer (1×) at a ratio of 10g/L, suspend it with a vortex shaker and add 0.5 mM PMSF.
  2. Lyse cells using ultrasound (3s pulse, 30 % amplitude, 8s pause).
  3. Centrifuge the crushing liquid for 10 minutes at 4°C, 12000rpm.
  4. Transfer the supernatant into a clean tube.
-

5. Aspirate the 20% ethanol in Ni-NTA agarose, and wash Ni-NTA agarose repeatedly with PBS buffer (1×) .
  6. For every 8 ml lysate, pipette 1 ml Ni-NTA agarose into a 10 ml tube.
  7. Place the test tube on a rotating incubator and incubate at 4°C for 90 minutes to make the protein and Ni-NTA agarose fully bound.
  8. Load the mixture onto the column and remove the column's bottom cap and collect flow-through in a 50 ml tube, which is placed in ice, until the mixture is completely loaded.
  9. Wash with 200ml Washing buffer A (table 1) and collect the flow-through with a 1.5ml EP tube (only one tube is sufficient), and mark the tube as "10mM".
  10. Wash with 200ml Washing buffer B (table 2) and collect the flow-through with a 1.5ml tube (only one tube is sufficient), and mark the tube as "50mM".
  11. Elute the protein with 1ml elution buffer (table 3) and collect the flow-through with a 1.5ml tube.
  12. Repeat the above steps 10 times and mark the tubes as "250-1 to 250-10".
  13. Wash with Washing buffer C (table 4), containing 500mM imidazole can clean the Ni-NTA agarose completely.
  14. Finally, wash with ddH<sub>2</sub>O and 20% ethanol successively, and store the Ni-NTA agarose in 20% ethanol.
  15. Analyze the collected flow-through sequentially in a SDS-PAGE.
  16. Load the purified proteins in a dialysis membrane and add an appropriate amount of rTEV protease (for Lba) or SUMO protease (for Lwa, Psm and Cca) to cleave off the His-MBP tag or the His-SUMO tag respectively.
-

17. Dialyse the sample in a large volume of PB Buffer (table 5) overnight at 4°C.
  18. Get rid of the cleaved off His-SUMO tag by repeating Ni-NTA agarose purification. This time however, the desired protein does not bind to the Ni-NTA agarose and can be collected in the first flow-through. Also elute the bound proteins to clean the Ni-NTA agarose for further use and also for SDS-PAGE analysis.
  19. Analyze the collected flow-through sequentially in a SDS-PAGE.
  20. Combine all the flow-through fluids containing Cas13 protein. Up-concentrate the sample to 3ml with a centrifugal filter (MWCO: 10 kDa) at 4500 rcf.
  21. Check purity in a SDS-PAGE.
  22. Measure the concentration via BCA Protein Quantitation Assay Kit.
  23. Store the protein at -80°C.
-