

Buffers and Media

Buffers

Citrate-Phosphate Buffer (McIlvaine-Buffer)¹

For 20 ml of buffer, mix 0.2 M Na₂HPO₄ and 0.1 M citric acid according to table:

Desired pH	0.2 M Na ₂ HPO ₄	0.1 M citric acid
2.2	0.40	19.60
2.4	1.24	18.76
2.6	2.18	17.82
2.8	3.17	16.83
3.0	4.11	15.89
3.2	4.94	15.06
3.4	5.70	14.30
3.6	6.44	13.56
3.8	7.10	12.90
4.0	7.71	12.29
4.2	8.28	11.72
4.4	8.82	11.18
4.6	9.35	10.65
4.8	9.86	10.14
5.0	10.30	9.70
5.2	10.72	9.28
5.4	11.15	8.85
5.6	11.60	8.40
5.8	12.09	7.91
6.0	12.63	7.37
6.2	13.22	6.78
6.4	13.85	6.15
6.6	14.55	5.45
6.8	15.45	4.55
7.0	16.47	3.53
7.2	17.39	2.61
7.4	18.17	1.83
7.6	18.73	1.27
7.8	19.15	0.85
8.0	19.45	0.55

2xSDS Sample Buffer

Ingredient	Volume in ml
1 M Tris pH 6.8	0.6
Glycerol	5
20 % SDS	1
Bromphenol blue	Small tip of spatula
1 M DTT	1
Water	2.4
Final	10

SDS Sucrose Buffer

Ingredient	Amount
20 % SDS	25 ml
Sucrose	15 g
Water	To 50 ml
Final	50 ml

DTT-Carbonat

Ingredient	Volume in μ l
1 M DTT	100
1 M Sodium-Carbonat	100
Water	800
Final	1000

Lämmli Running Buffer

Ingredient	Volume in ml
250 mM Tris	100
1.34 M Glycine	100
20 % SDS	5
0.5 M EDTA	2
Water	To 1 l
Final	1 l

Media

Revised TAP Medium

See Kropat *et. al*, 2011²

TAP-Agar Plates

Add 15 g of select agar to 1 l of revised TAP medium. Autoclave for 1 h.

LB-Medium

Add 10 g Tryptone, 5 g yeast extract and 10 g sodium chloride together. Fill up to 1 l with aqua dest.

Autoclave for 1 h.