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Purpose

To compare the half-life of Green fluorescent protein (GFP) modified by TJU with that of unmodified GFP.

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Material & Equipment

Microplate tester, ultrasonic splitting tester, sterilizing pot, liquid culture medium, Sterile head, pipette gun, ultra clean table. *Saccharomyces cerevisiae* strains ① that have modified GFP integrated on their genomes and can express CRISPR/Cas9 which can complete genome cutting purposes. *Saccharomyces cerevisiae* strains ② that have GFP integrated on their genomes and can express CRISPR/Cas9 with the same function of strains ①.

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Procedure

A. Detection of fluorescence intensity by microplateReader

Inoculate ① and ② respectively in liquid medium, and set the speed and temperature as the optimal conditions of *S. cerevisiae*.

At the incubation time of 0/12/24/48/72h, 1ml sample was taken in a clean bench with sterile head (strictly sterile operation, avoid bacterial contamination), and 200 μ l sample was added to the 96-well plate after sampling. Excitation light was set for the microplate reader program, and absorbance was measured at appropriate emission wavelength and recorded.

B. SDS-PAGE

Inoculate ① and ② respectively in liquid medium of at least 200ml, and set the speed and temperature as the optimal conditions of *S. cerevisiae*.

The precipitate was recovered at 4000rpm, centrifuged at 4℃ for 10min.

Add an appropriate amount of buffer solution (10% glycerol, 20Mm HEPES, 500mM NaCl, adjust the pH value to the optimal conditions of GFP, generally 7.5) to the precipitation.

Ultrasonic schizotometer lyses thallus. Set to 600W, ultrasound 5s, 5s interval, a total of 40min.

Centrifuging at 15000rpm, 4℃ for 20min(Repeat 3 times). Take the supernatant and reserve.

SDS-PAGE samples were prepared with 20 μ l supernatant + 5 μ l Loading Buffer and heated at 98℃ for 15min to complete denaturation.

Sample loading, swim for 45min at 150V.

After electrophoresis, stain with staining solution (Coomassie bright blue R250 1.2g/L; Methanol: water: acetic acid (5:4:1)) at 60℃ for 8min.

After dyeing, decolorization solution (methanol: water: acetic acid (5:4:1)) was used to decolorize until a clear band was visible.

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Data processing

The fluorescence intensity was measured with a microplate reader and the data were derived to draw a histogram

SDS-PAGE can be qualitative, only semi-quantitative. When the protein concentration is high enough, fluorescence intensity can be measured using a developer.

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Precautions

The above conditions need to be adjusted and optimized according to the properties of *S. cerevisiae*.

The excitation light and emission light for fluorescence detection should be adjusted according to the actual situation of GFP used.

Do not use manual dyeing method if there is a glue dyeing machine.