

Gel Electrophoresis

Purpose

Agarose gel electrophoresis will be used in order to determine the success of PCR reactions, restriction digestions, ligations and mini-prepped plasmids after cloning and transformation.

Materials

- ❖ DNA (Label with name)
- ❖ Nuclease-free H₂O
- ❖ 4 µL Loading Dye (6X)

- ❖ PCR tubes
- ❖ 250 mL Erlenmeyer flask
- ❖ DNA samples for loading onto gel (PCR, Restriction digestion product, Ligation product, Mini-prepped plasmid)
- ❖ Nuclease-free H₂O
- ❖ Double-distilled H₂O
- ❖ RedSafe Nucleic Acid Staining Solution (20,000X)
- ❖ 1X TAE Buffer - 40 mM Tris, 20 mM acetic acid, 1 mM Na₂EDTA
- ❖ Agarose
- ❖ Gel Loading Dye, Purple (6X)
- ❖ 1 kb Plus DNA Ladder
- ❖ 50 bp DNA Ladder

Equipment

- ❖ Top-loading balance
- ❖ Bio-Rad Gel imaging station
- ❖ Electrophoresis power source
- ❖ Agarose gel casting apparatus
- ❖ Microwave Oven

Procedure

Preparing the Agarose Gel

- Weigh the appropriate amount of agarose necessary based on the agarose concentration determined above. For a 1.5% gel, weigh 1.5 g of agarose.
- Add 100 mL of 1X TAE to a 250 mL Erlenmeyer. Add the agarose powder. Swirl to mix.

- Place the Erlenmeyer inside the microwave oven. Microwave on high for 2 minutes. Watch the contents of the microwave closely to ensure that the contents do not boil over.
 - a. If the solution is not complete clear and free of agarose particles, then microwave on high for another 30 seconds.
- Place the hot Erlenmeyer flask on a heat-resistant surface and allow it to cool.
- When the surface of the Erlenmeyer is warm, but not hot to the touch, add 5 μ L of RedSafe Nucleic Acid Staining Solution (20,000X) to the agarose solution. Swirl the flask gently to mix and avoid forming bubbles.
- Place an agarose gel tray inside the gel caster and tightly seal the caster against the open sides of the gel tray. Place a comb inside the casting apparatus at the appropriate grooves on the side of the tray.
 - a. If you poured water into the casting unit, no liquid should leak from any side.
- Pour the agarose solution into the gel tray. You should be immersing between half and a quarter of the comb length in agarose.
- Leave on a heat-resistant surface and allow to cool and solidify. This will take approximately 20 minutes.
- When the gel is cool, carefully pull out the comb.
- Place the solidified gel kept in the casting tray inside the electrophoresis chamber. The wells should be located at the BLACK negative end of the apparatus.
- Pour 1X TAE buffer on top of the gel into the electrophoresis chamber so that it covers the gel with about 3 mm of liquid.

Sample Preparation, Loading and Gel Electrophoresis

- Label each PCR tube with the name of sample being tested.
- For each tube, add the volumes of nuclease-free water, DNA sample and 6X loading dye calculated above. Pipette up and down gently to mix.
- Load the first well of the gel with 5 μ L of the ladder chosen above.
- Load the wells of the gel with the samples prepared in PCR tubes.
 - a. Do not move too slowly when loading the wells, otherwise the first samples will begin to diffuse as the last samples are being loaded.
- Place the lid on the gel electrophoresis chamber, taking care that the red electrode matches the red end on the chamber.
- Turn on the electrophoresis power source. Plug the electrodes connected to the gel electrophoresis chamber lid to the electrophoresis power source. Please consult SOP #17 Electrophoresis power supply for additional details prior to usage.
- Run the gel electrophoresis at a constant voltage of 120 V for 60 minutes. Monitor the progression of the gel to ensure that the band representing the loading dye does not run off the gel.
 - a. Stop the electrophoresis if this band is about to run off.
 - b. Extend the duration of electrophoresis if this band does not pass ~75% of the gel.
- Stop the run and turn off the gel electrophoresis power supply.

Gel Imaging

- Remove the gel tray containing the gel from the gel electrophoresis unit.
- Open the door of the transilluminator (BioRad Gel Doc XR+) and place the gel on the transilluminator surface. Close door.
 - Consult SOP #19 Bio-Rad Gel Imager for more information
- Select the “GelRed” protocol and follow computer instructions for changing the filter settings. Run the protocol and either utilize automatic exposure or manual exposure.
- Acquire the image and save as jpeg in team folder. Send image through email on the laptop.

Safety Precautions

- Heat, electrical current, UV light
- Always wear safety goggles, lab coat and disposable gloves.
- Wash hands in the designated hand-washing sink once finished.
- Wear appropriate skin and eye protection when working with UV light.
- Ensure all switches and indicators are in proper working condition.
- Ensure power cords and leads are undamaged and properly insulated.
- Wear insulated work gloves when heating and transferring agarose gel. Do not cap the agarose gel container when heating.

Reference

- <https://www.thermofisher.com/content/dam/LifeTech/global/brands/Documents/1114/general-recommendations-dna-electrophoresis.pdf>
- IBEHS 2P03 Laboratory Manual
- SOPs supplied by ETB 311