

Gibson Assembly

Gibson Assembly

Introduction

Gibson Assembly facilitates plasmid cloning in a one-pot reaction. In this context, a single fragment has been joined to the backbone plasmid.

Materials


- › Backbone plasmid
- › DNA Insert
- › Gibson Assembly Master Mix (2X)
- › Deionized H₂O

Procedure

Gibson Assembly

1. Prior to reaction assembly, DNA samples must be diluted to a determined concentration that contemplates 0.02–0.5 pmols of total DNA quantity.
2. Mix the previously noted components as follows:
 - x µL DNA insert
 - x µL Backbone Plasmid
 - 10 µL Gibson Assembly Master Mix
 - 10-2x µL Deionized H₂O (or up to 20 µL)
3. Incubate samples in a thermocycler at **50°C for 15 minutes**
4. Store samples on ice or at -20°C for subsequent transformation

Further detail of the steps of Gibson Assembly cloning [2]



Bibliography

- [1] New England Biolabs. (2012, 25 septiembre). *Gibson Assembly® Master Mix – Assembly (E2611)* | NEB. <https://international.neb.com/protocols/2012/09/25/gibson-assembly-master-mix-assembly>
- [2] Cutts, E., & Vannini, A. (2018). Troubleshooting biGBac: a practical guide v1. *protocols.io*. Published. <https://doi.org/10.17504/protocols.io.q3ydpw>