

In Vitro gRNA- Cas12a Digestion Protocol

In Vitro gRNA-Cas12a Digestion

Introduction

Protocol obtained from NEB, specific for the purified protein that was ordered.

EnGen Lba Cas12a (Cpf1) from Lachnospiraceae bacterium ND2006 is a site-specific DNA endonuclease guided by a single 41-44 nucleotide guide RNA (gRNA) (1). Targeting requires a gRNA complementary to the target site as well as a 5' TTTN protospacer adjacent motif (PAM) on the DNA strand opposite the target sequence. Cleavage by EnGen® Lba Cas12a occurs ~18 bases 3' of the PAM and leaves 5 nucleotide 5' overhanging ends. EnGen Lba Cas12a has Simian virus 40 (SV40) T antigen nuclear localization sequences (NLS) at both the N and C-termini of the protein.

Materials

- › EnGen Lba Cas12a (Cpf1) (NEB #M0653)
- › 10X NEBuffer r2.1 Reaction Buffer
- › Nuclease-free water
- › Proteinase K, Molecular Biology Grade (NEB #P8107)
 - › Subtilisin-related serine protease hydrolyzes a variety of **peptide bonds** and is frequently used to cleanup enzymatic reactions or cell lysates. Probably won't be needed if we amplify the DNA target beforehand.
- › Guide RNA containing the targeting sequence in the region of interest --> our construct
 - › AmpR
 - › ChloR
 - › EryR
 - › KanR
 - › SpecR
- › DNA substrate containing the target sequence --> lab plasmids with AR sequences
 - › The substrate DNA can be circular or linearized plasmid, PCR products, or synthesized oligonucleotides
- › Optional Materials
 - › Apparatus and reagents for DNA fragment analysis --> Agarose gel electrophoresis apparatus

Procedure

Before starting

1. Wear gloves and use **nuclease-free tubes** and reagents to avoid RNase contamination.
2. Reactions are typically **30 µl** but can be scaled up as needed. Reactions should be assembled in nuclease-free microfuge tubes or PCR strip tubes.
3. It is essential to keep the **molar ratio** of Cas12a and gRNA per target site at **10:10:1** or **higher** to obtain the best cleavage efficiency. A calculator can be found here: <https://nebiocalculator.neb.com/#!/ligation>
4. Prepare **300 nM gRNA** by diluting the stock with nuclease-free water on ice.

5. Prepare **30 nM substrate DNA** with a single target sequence by diluting the stock with nuclease-free water on ice.
6. If planning to use higher concentration EnGen Lba Cas12a (NEB #M0653T) for *in vitro* digestion of DNA, the enzyme can be diluted to **1 μ M** in EnGen Lba Cas12a Diluent prior to the reaction assembly.

In vitro Assay

7. Assemble the reaction at **room temperature** in the following order:

	A	B
1	COMPONENT	AMOUNT
2	Nuclease-free water	20 μ l
3	NEBuffer r2.1 Reaction Buffer (10X)	3 μ l
4	300 nM gRNA	3 μ l (30 nM final)
5	1 μ M EnGen Lba Cas12a (Cpf1)	1 μ l (~30 nM final)
6	Total Reaction Volume	27 μl

8. **Pre-incubate** for 10 minutes at 25°C.

9. **Add** 3 μ l of 30 nM substrate DNA (30 μ l final volume).

10. **Mix** thoroughly and pulse-spin in a microfuge.

11. **Incubate** at 37°C for 10 minutes.

12. NOT: Add 1 μ l of Proteinase K (NEB #P8107) to each sample, Mix thoroughly, and pulse-spin in a microfuge. (Stop reaction, Cas12a digestion)

13. **Incubate** at room temperature for 10 minutes.

14. **Analyse** with the detection protocol [2].

Bibliography

1. New England Biolabs. (n.d.-b). In vitro digestion of DNA with EnGen® Lba Cas12a (Cpf1) (M0653). Retrieved October 17, 2021, from Neb.com website:
<https://international.neb.com/protocols/2017/12/19/in-vitro-digestion-of-dna-with-engen-lba-cas12a-cpf1-neb-m0653>
2. 1. (N.d.). Retrieved October 17, 2021, from Windows.net website:
https://sfvideo.blob.core.windows.net/sitefinity/docs/default-source/user-guide-manual/dnasealert-detection-manual.pdf?sfvrsn=34133407_8