

5. PCR & Agarose Gel Electrophoresis – Colony PCR

Material

Green Taq Mix
Primer
Double Distillation Water (ddH₂O)
Eight tube centrifuge
EP Tube
Tweezers
plastic wrap
Shaker
PCR Eight-Tube
DNA Marker

Step(Take 5 tubes of eight tube centrifuge for example)

- 1) Pick bacteria
 - ① Add 1mL Amp LB to 1.5mL EP Tube.
 - ② Use a small tip with tweezers to pick out a single clone and put them into the EP Tube.
 - ③ Put the EP tube on the Rank and wrap it with plastic wrap.
 - ④ Insert it into the shaker obliquely, shake for 2h.
- 2) Colony PCR (after shaking for 2h)
 - ① Adding 6.6 Primer Mix(5 μ l), 41.3 μ l 2 x Green Taq, 29.1 μ l of ddH₂O to the first tube of eight-tube, then draw 14ul from the first tube to the next 4 tubes.
 - ② Centrifuge the eight-tube in a centrifuge for several seconds.
 - ③ Take the eight-tube into Clean bench.
 - ④ Add 1 μ l Bacteria Liquid.
 - ⑤ Set the protocol as follow: begin at 95°C for 3 mins, then keep 95°C for 15 secs for denaturation, decrease to 55°C for 15 secs, 72°C for 1 min and repeat that cycle 30 times, 72°C for 1 min, finally maintain 4°C infinity.
 - ⑥ Adding all the samples to the hole. Run the gel at 120V for about 25 minutes. Check the result under the Blue Light Gel Imager.
 - ⑦ Select the right System for Shaking bacteria and Plasmid Extraction.

Note

- ① The lid of the PCR Eight-tube is placed in the ultra-clean workbench to make sure we add Bacteria Liquid at the ultra-clean workbench.
- ② As we use Green Taq Mix, we needn't add Loading when Agarose Gel Electrophoresis.