

***Synechocystis* sp. PCC 6803 Transformation Protocol**

Extracted from: Protocol based on protocol from the University of Turku received from Dr. Julie Zedler. Version: September 2017

Materials:

- *Synechocystis* target strain culture (healthy, logarithmic phase)
- 5 µg of purified plasmid per transformation
- fresh BG-11 medium
- BG-11 (NO antibiotics) plates
- BG-11 + 20% Antibiotics plates (depend on what the plasmid is you are using)

Protocol:

1. Day 1
 - a. Pellet 5 - 10 mL of a healthy target strain (logarithmic phase, OD₇₅₀ around 1-3) per transformation (4000 rpm, 10 min, RT in falcon tubes) (don't forget a negative control!!!).
 - b. Discard supernatant and resuspend cells in 100 µL fresh BG-11 per 5-10 mL culture used.
 - c. Transfer 100 µL of cells into clean eppy, add approximately 5 µL (1 µg) of plasmid DNA. Mix the suspension by pipetting carefully.
 - d. Incubate the mixture shaking O/N at 30°C.
 - e. Remark: I usually put the eppys in a brown 5 mL falcon tube and put that in the shaking incubator.
2. Day 2
 - a. Plate the cell/DNA suspension straight on BG-11 plates (no antibiotics!).
 - b. Spread cells and leave plates to dry properly in the sterile hood.
 - c. Incubate O/N in the 30°C incubator with light protection (cover with some tissue).
3. Day 3
 - a. Remove light protection and incubate in full light until the next day.
4. Day 4
 - a. Wash cells carefully off the plate with approx. 400 µL fresh BG-11.
 - b. Remark: I find it best to hold the plate vertical with one hand and repeatedly flush cells off to the bottom by pipetting.
 - c. Transfer suspension on BG-11 plates with 20% antibiotics.

- d. Spread cells and leave plates to dry properly in the sterile hood.
 - e. Incubate in full light at 30° C for approx. 7 to 14 days.
5. After colonies are visible....
- a. Re-streak single colonies on full concentration antibiotics BG-11 plates.
 - b. Check for example by PCR. (Full segregation can take quite a while, keep on antibiotics).