

PCR Using Q5 High-Fidelity DNA Polymerase

Prepare MasterMix:

Per PCR-Tube:

- 5 μ L 5x Q5 Reaction Buffer
- 0,5 μ L 10 mM dNTPs
- 0,25 μ L High-Fidelity Q5 Polymerase
- 15,75 μ L MQ (to get to 25 μ L with DNA and Primers)

-> snip then spin down

-> transfer 21,5 μ L to each PCR-tube

Add DNA and Primers:

Per PCR tube:

- 1 μ L Template DNA
- 1,25 μ L 10 μ M forward Primer
- 1,25 μ L 10 μ M reverse Primer

PCR-Program:

1. Initial Denaturation : 30 s at 98°C
2. Annealing and Extension: -> 10 s at 98°C
(25-35 cycles) -> 30 s at 50-72°C*
-> 20-30 s at 72°C per 1000 bp
3. Final Extension : 30 s at 72°C
4. Hold : 8°C

* use [NEB Tm Calculator](#)