

# Sequential Double Digest

---

## Introduction

This is the Sequential Double Digest Protocol with Standard Restriction Enzymes. If there is no buffer in which the two enzymes exhibit > 50% activity, this sequential digest can be performed.

More information from NEB can be found [here](#).

Double Digests can be designed using [NEB's Double Digest Finder](#).

[NEBcloner](#) will help guide your reaction buffer selection when setting up double digests.

See the [NEBuffer Activity/Performance Chart with Restriction Enzymes](#) for the incubation temperatures.

## Materials

- › DNA **1  $\mu$ g**
- › [NEBuffer](#)
  - › 1X
- › [NEB Restriction Enzymes](#)
- › Deionized Water

## Procedure

### Sequential Double Digest

- ✓ 1. Set up the following reaction using the restriction endonuclease that has the lowest salt concentration in its recommended buffer (total reaction volume **50  $\mu$ l**).
- ✓ 2. Set up the following digest reaction on ice

Table1				^
	A	B	C	
1		Volume (μl)		
2	Buffer (10x)	5		
3	DNA *	Input Volume for ng		
4	Restriction Enzyme #1 **	1		
5	Deionized Water (μl)	1		
6	Total Volume (μl)	#VALUE!		

\*A 50 μl reaction volume is recommended for digestion of 1 μg of substrate.

\*\* Restriction Enzyme, 10 units is sufficient, generally **1 μl** is used

\*\*\*The enzyme should be the last component added to reaction

- ✓ 3. Mix components by pipetting the reaction mixture up and down, or by "flicking" the reaction tube.
- ✓ 4. Quick ("touch") spin-down in a microcentrifuge. Do not vortex the reaction.
- ✓ 5. Incubate for 1 hour at the enzyme-specific appropriate temperature.

01:00:00



Can be decreased to 5-15 minutes by using a [Time-Saver™ Qualified Restriction Enzyme](#)

See the [NEBuffer Activity/Performance Chart with Restriction Enzymes](#) for the incubation temperatures.

- ✓ 6. Adjust the salt concentration of the reaction (using a small volume of a concentrated salt solution) to approximate the reaction conditions of the second restriction endonuclease.&nbsp;

Alternatively, a spin column can be used to isolate the DNA prior to the second reaction.

- ✓ 7. Add the second enzyme.

\*\* Restriction Enzyme #2, 10 units is sufficient, generally **1 μl** is used

- ✓ 8. Mix components by pipetting the reaction mixture up and down, or by "flicking" the reaction tube.
- ✓ 9. Quick ("touch") spin-down in a microcentrifuge.&nbsp;Do not vortex the reaction.
- ✓ 10. Incubate for 1 hour at the enzyme-specific appropriate temperature.

Can be decreased to 5-15 minutes by using a [Time-Saver™ Qualified Restriction Enzyme](#)

See the [NEBuffer Activity/Performance Chart with Restriction Enzymes](#) for the incubation temperatures.