

Ligation Protocol WITH T4 DNA Ligase (M0202)

Introduction

Please see the [NEB website](#) for supporting information on this protocol.

Materials

- › 10X T4 DNA Ligase Reaction Buffer
- › T4 DNA Ligase
- › Vector DNA (4kb)
- › Insert DNA (1kb)
- › Nuclease-free water

Procedure

Set up the T4 DNA Ligase Reaction

Note: T4 DNA Ligase should be added last. The table shows a ligation using a molar ratio of 1:3 vector to insert for the indicated DNA sizes. Use [NEB calculator](#) to calculate molar ratios.

- ✓ 1. Thaw the T4 DNA Ligase Buffer and resuspend at room temperature.

Tip: Aliquote the 10x buffer less concentrated so when thawing, the DTT gets soluble more easily.

- ✓ 2. Set up the following reaction in a microcentrifuge tube on ice:

Table1			^
	A	B	
1	Component	Volume (µl)	
2	10X T4 DNA Ligase Buffer	2	
3	Vector DNA: 50 ng (0.020 pmol)		
4	Insert DNA: 37.5 ng (0.060 pmol)		
5	Nuclease-free water	17	
6	T4 DNA Ligase	1	
7	Total	20	

- ✓ 3. Gently mix the reaction by pipetting up and down and microfuge briefly.
- ✓ 4. For cohesive (sticky) ends, incubate at 16°C overnight or room temperature for 10 minutes. For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours.

- ✓ 5. Heat inactivate at 65 degrees C for 10 minutes.

00:10:00



- ✓ 6. Chill on ice and transform 1-5 μ l of the reaction into 50 μ l competent cells.

Use 25 μ L DH5 α cells, and add 2 μ L of reaction mixture.