

EXPAR - DNA Protocol (65oC)

Introduction

An isothermal amplification method for microRNAs.

End Product: DNA

Reaction Temperature: 65°C

The protocol is based on Ellie Mok *et al.*; [Comprehensive evaluation of molecular enhancers of the isothermal exponential amplification reaction](#).

Materials

› Materials

- › Template DNA (instead of primers) for each microRNA
- › microRNAs
- › Nb.BtsI (nicking enzyme)
- › rCutSmart™ Buffer (buffer for Nb.BtsI)
- › Bst 2.0 DNA Polymerase
- › Isothermal Amplification Buffer (buffer for Bst 2.0 DNA Polymerase)
- › Deoxynucleotide (dNTP) Solution Mix
- › SYBR Green I
- › PCR H₂O

› Equipment

- › Pipettes & tips
- › Eppendorfs
- › Centrifuge
- › Incubator
- › Real Time PCR machine

›

Procedure

1. Prepare Mix A (for each microRNA):

Mix A			
	A	B	C
1	Material	Final Concentration	Volume
2	PCR H ₂ O		up to 20 ul final volume (for both Mix A & Mix B)
3	Template DNA	150 nM	0.5 ul
4	microRNA	10 nM	0.5 ul
5	dNTPs	250 mM	0.5 ul
6	rCutSmart Buffer	0.5x	1 ul

2. Prepare mix B (for each microRNA):

Mix B			
	A	B	C
1	Material	Final Concentration	Volume
2	Isothermal Amplification Buffer	1x	2 ul
3	Bst 2.0 DNA Polymerase	0.4 U/ul	1 ul
4	Nb. BtsI	0.5 U/ul	1 ul
5	SYBR Green	0.5x	1 ul

3. Incubate mix A to 56°C for 10 minutes (hybridisation of template and microRNA)

4. Add mix B to mix A while on ice

5. Prepare the necessary negative controls by following the above-mentioned procedure but without adding microRNA to Mix A and by adjusting the right volume of water

6. Centrifugation of the samples for 3 sec. at max speed

7. Load 20 ul of the samples to the wells of the q-PCR plate while on ice

8. Put the plate on the q-PCR machine and run the designed protocol (65°C for 20 minutes)