

# Promoters Optimization

05/09/2021

**Cloning of optimized promoters for *E. coli* and *B. subtilis* into the original plasmid, by replacing the original promoter with PCR amplification and Gibson Assembly:**

PCR amplification of mCherry plasmid and exclusion of the p43 promoter:

First, six and five optimized promoters for *B. subtilis* and *E. coli*, respectively, were chosen and labeled according to the table below:

	A	B	C	D	E
1	<u>Bacteria</u>	<u>Code</u>	<u>Name by modeling team</u>	<u>Primers' names</u>	<u>gBlock</u>
2	B.subtilis	A1-1	ribosomal protein s10	rpsj Fwd (Rev same as rpsj syn!)	IDT
3	B.subtilis	A1-2	ribosomal protein s10 synthetic version	rpsj syn Fwd + Rev	IDT
4	B.subtilis	D1	L-lactate permease	lactate Fwd + Rev	twist (D1)
5	B.subtilis	A2	ribosomal protein L4	L4 Fwd + Rev	twist (A2)
6	B.subtilis	B2	rpld2 (synthetic)	rpld syn Fwd (Rev same as A2 [L4])	twist (B2)
7	B.subtilis	E1	lctp4	lctp syn (primers same as D1 [lactate])	twist (E1)
8					
9	E.coli	C2	cAMP phospho	cAMP Fwd + Rev	twist (C2)
10	E.coli	D2	cpda synthetic	cpda syn (same primers as C2 [cAMP])	twist (D2)
11	E.coli	M1	RNA polymerase transcription factor DK (dksa)	dksa Fwd + Rev	IDT
12	E.coli	T1	16S rRNA pseudo	rsua Fwd + Rev	IDT
13	E.coli	K1	rpma	rpma Frw+Rev	IDT

The mCherry vector was PCR-amplified with corresponding primers containing overlapping tails for its gBlock (=promoters),(see Primers document), in the following thermocycling conditions:

- Initial Denaturation: 98°C, 3 sec
- 30 cycles:
  1. 98°C, 10 sec
  2. 57°C, 30 sec
  3. 72°C, 3:50 min
- Final extension: 72°C, 2min
- Hold: 10°C
- Total volume: 50ul

10ul from each PCR product was examined by a DNA electrophoresis gel as below:  
The rest of the volume was kept at -20°C.



Loading scheme: T1, A1S, K1, A2, A1, B2, M1, C2, D1, E1 & D2.

09/09/2021

The rest of the PCR products were treated with DpnI in a total volume of 50ul, 10% CutSmart buffer for 1hr at 37°C.

Since PCR amplification for promoters cloning resulted in unspecific bands (see the gel above), we ran the total volume of PCR products following degradation by DpnI in a wide-wells agarose gel, and relevant bands (~7400bp) were extracted by GenElute Gel Extraction Kit from Sigma, according to manufacturer's instructions. Then, concentrations of linearized plasmids and gBlocks (promoters from IDT and Twist= inserts) were measured by Nanodrop:

Concentrations of PCR products and gBlockes and calculation for Gibson Assembly reaction volumes								
	A	B	C	D	E	F	G	H
1	PCR product	PCR Conc	Gibson	gBlock	Conc	Gibson	Vector+insert	UPW
2	A2	18.1	6.0	A2	8.7	2.5	8.5	1.5
3	B2	5.5	9.1	B2	11.6	0.5	9.6	0.4
4	C2	24.4	4.4	C2	11.1	2.0	6.4	3.6
5	D1	23	4.7	D1	12.3	1.8	6.5	3.5
6	D2	25.7	4.2	D2	11	2.0	6.2	3.8
7	E1	19.2	5.6	E1	10.4	2.1	7.7	2.3
8	K1	28.3	3.8	K1	15.8	1.4	5.2	4.8
9	M1	17.7	6.1	M1	13.3	1.7	7.8	2.2
10	T1	13.5	8.0	T1	15.3	1.4	9.4	0.6
11	A1	11.5	9.4	A1	19.7	1.1	10.5	0
12	A1S	14	7.7	A1S	14.7	1.5	9.2	0.8
≡								

The table above shows measured concentrations of vectors and inserts, as well as the amount in  $\mu\text{l}$  needed for Gibson assembly reaction, in a molar ratio of 1:2 in favor of the insert.

13/09/2021

#### Gibson Assembly for promoters:

Gibson assembly was performed as follows:





Vector (from 09.09.21): gBlock ratio was 1:2, in which 0.05pmols of gBlock and 0.025pmols of vector were added to 10ul of HiFi DNA Assembly Master Mix (NEB) and deionized water to make 20ul. Samples were incubated at 50°C for 1hr. Then, 3ul from each reaction was used for bacterial transformation according to NEB's protocols. 100ul of cells were spread over an LB- agar + chloramphenicol 10ug/ul plates and incubated ON at 37°C.

17/09/2021

#### Screening new primer-positive colonies by colony PCR:

Colonies were picked for colony PCR assessment in order to find positive colonies having the new insert. Since the size of the vector with the insert is larger by ~360bp than the original vector, colony screening was based on the vector's size instead of insert-specific primers.

Colony PCR was performed on four or five colonies from each transformation, according to the protocol. Additionally, bacterial suspensions from picked colonies were streaked over an LB-agar plate.





18/09/2021

Insert-positive colonies were transferred to ON Starters, except A1S and B2 that failed to grow.

19/09/2021

Plasmids were prepared from all ON starters with PureYield plasmid miniprep system (Promega). And plasmids concentrations were measured via NanoDrpo and examined in DNA electrophoresis gel:

Sample	conc. (ng/ul)	260/280	270/230
M-1	435.8	1.88	2.3
M-4	438.5	1.89	2.28
M-5	472	1.89	2.34
A1-2	226.8	1.88	2.25
A1-4	176.3	1.87	2.19
A1-5	262.2	1.89	2.28
TAI-D1	238.8	1.83	1.68
TAI-D2	223.4	1.88	2.27
TAI-D3	239	1.88	2.25
T-1	433.8	1.89	2.29
T-3	378	1.87	2.19
T-4	442.3	1.87	2.14
E-1	474.9	1.88	2.29
E-2	412.7	1.87	2.27
E-3	600.8	1.91	2.32
K-1	630.3	1.9	2.3
K-2	372.3	1.87	2.3
K-4	263.5	1.88	2.15
D1-2	390.8	1.891	2.33
D1-3	375.2	1.88	2.27
D1-4	251.1	1.84	2.01
C-1	382.6	1.88	2.25
C-4	334.7	1.86	2.11
C-5	380.6	1.88	2.22
A2-1	453.3	1.9	2.3
A2-4	373.8	1.87	2.07
A2-5	364.1	1.88	2.2
D2-2	367.3	1.87	2.19
D2-3	371.8	1.87	2.06
D2-5	232.1	1.87	2.11



Order in gel from left to right is identical to the list in the table.

22/09/2021

Plasmids were sent for sequencing in hylabs.

23/09/2021

From each plasmid preparation, at least one was positive for the new promoter version. One plasmid from each version was transformed into *E. coli* and *B. subtilis* following the cloning protocols.

24/09/2021

One colony from each transformation and bacterial species was picked for ON starter.

25/09/2021

Glycerol stock was prepared from each starter.

03/10/2021

#### Assay for promoters:

Bacteria containing the new promoter versions (from glycerol stocks) were streak over an LB-agar plate with 10mg/ml chloramphenicol.

04/10/2021

Bacteria containing the new promoter versions (from glycerol stocks) were streak over an LB-agar plate with 10mg/ml chloramphenicol, for a second biological repeat.

Colonies picked to ON starters from plates of the previous day.

05/10/2021

The assay was performed according to the protocol. Briefly, bacteria from ON starters were washed with PBSx1, and LB was replaced with M9 medium. Following 4hrs incubation, OD600nm was measured and bacteria were diluted to OD600nm of 0.2, according to the table below:

Bacteria	Sample	D600nm	Bacterial suspension for dilution	M9 for dilution
<i>E. coli</i>	ECE	0.817	171	529
	mCherry	1.149	122	578
	A2	0.829	169	531
	C	0.752	186	514
	D1	1.059	132	568
	D2	1.003	140	560
	E	1.035	135	565
	K	1.054	133	567
	M	0.794	176	524
	T	1.279	109	591
<i>B. Subtilis</i>	ECE	0.805	174	526
	mCherry	0.731	192	508
	A2	0.657	213	487
	C	0.747	187	513
	D1	0.366	383	317
	D2	0.925	151	549
	E	0.273	513	187
	K	0.51	275	425

	<b>M</b>	0.622	225	475
	<b>T</b>	1.004	139	561

Then, 200 µl of each bacterial sample were loaded into 96 well plates in triplicates, as described in the plating scheme:

Well6												
	1	2	3	4	5	6	7	8	9	10	11	12
A	Coli-ECE	Coli-ECE	Coli-ECE	Coli-M	Coli-M	Coli-M	Bac-ECE	Bac-ECE	Bac-ECE	Bac-M	Bac-M	Bac-M
B	Coli-mCherry	Coli-mCherry	Coli-mCherry	Coli-T	Coli-T	Coli-T	Bac-mCherry	Bac-mCherry	Bac-mCherry	Bac-T	Bac-T	Bac-T
C	Coli-A2	Coli-A2	Coli-A2	Blank	Blank	Blank	Bac-A2	Bac-A2	Bac-A2			
D	Coli-C	Coli-C	Coli-C				Bac-C	Bac-C	Bac-C			
E	Coli-D1	Coli-D1	Coli-D1				Bac-D1	Bac-D1	Bac-D1			
F	Coli-D2	Coli-D2	Coli-D2				Bac-D2	Bac-D2	Bac-D2			
G	Coli-E	Coli-E	Coli-E				Bac-E	Bac-E	Bac-E			
H	Coli-K	Coli-K	Coli-K				Bac-K	Bac-K	Bac-K			

Moreover, colonies picked to ON starters from plates of the previous day.

06/10/2021

The assay was performed as mention above.

OD600nm measurements and dilution to OD600nm=0.2 calculations:

<b>Bacteria</b>	<b>Sample</b>	<b>D600nm</b>	<b>Bacterial suspension for dilution</b>	<b>M9 for dilution</b>
<b><i>E. coli</i></b>	<b>ECE</b>	0.893	157	543
	<b>mCherry</b>	1.676	84	616
	<b>A2</b>	0.889	157	543
	<b>C</b>	1.188	118	582
	<b>D1</b>	0.949	148	552
	<b>D2</b>	1.133	124	576
	<b>E</b>	1.371	102	598
	<b>K</b>	1.055	133	567
	<b>M</b>	1.808	77	623
	<b>T</b>	1.046	134	566
<b><i>B. Subtilis</i></b>	<b>ECE</b>	0.866	162	538
	<b>mCherry</b>	0.967	145	555
	<b>A2</b>	1.086	129	571
	<b>C</b>	0.749	187	513
	<b>D1</b>	0.889	157	543
	<b>D2</b>	0.86	163	537
	<b>E</b>	1.153	121	579
	<b>K</b>	0.802	175	525
	<b>M</b>	0.894	157	543
	<b>T</b>	0.814	172	528

Then, 200 µl of each bacterial sample was loaded into 96 well plates in triplicates, as shown in the plating scheme above.