

# Cloning using restriction enzymes

## Protocoles

Fragment amplification via PCR.....	p2
Insert digestion.....	p4
Plasmid digestion.....	p6
Ligation.....	p7
BL21(DE3) E. coli heat shock transformation.....	p8



# Fragment amplification via PCR

## Introduction

This is the protocol for fragment amplification with NEB® Hot Start Taq 2X Master Mix to get your insert ready for molecular cloning using restriction enzymes..

## Material

- Reverse and Forward primers
- Template DNA (100-200ng/μL)
- NEB® Hot Start Taq 2X master mix
- Nuclease-free water
- Bucket of ice
- Promega Wizard® SV Gel and PCR Clean-Up System kit

## Procedure

1- Add to a sterile thin-walled PCR tube on ice:

component	50 μL reaction
20 μM Forward Primer	1.5μL
20 μM Reverse Primer	1.5μL
Template DNA	2μL
NEB Hot Start Taq 2X master mix	25μL
Nuclease-free H2O	20μL

Gently mix the reaction and perform a quick spin to collect all liquid to the bottom of the tube.

2- Transfer PCR tubes to a PCR machine and run the following thermocycling protocol:

Step	Temperature (°C)	time
Initial denaturation	95°C	30 secondes
35 cycles	95 °C 55°C 68 °C	30 secondes 30 secondes 50 secondes
Final extension	68 °C	5 minutes



Hold	10°C	infinite
------	------	----------

3- Perform a PCR clean-up using the Promega Wizard® SV Gel and PCR Clean-Up System.

Add 50µL of Membrane Binding Solution to your PCR amplification product and follow the 5.A. DNA Purification by Centrifugation part of the protocol given in the kit.

4- Nanodrop your sample to make sure of the purity of your sample and to know the concentration. Your insert can now be digested.



# Insert digestion

## Introduction

This is the protocol for fragment digestion to get ready for molecular cloning with restriction enzymes.

## Material

- Ice bucket
- DNA insert
- NEB® rCutSmart Buffer
- NEB® 3.1 Buffer
- NEB® XhoI
- NEB® XhoI
- NEB® XhoI
- NEB® EcoRI
- Nuclease-free water
- TAE 1X
- Low Melting Agarose powder
- SYBR™ Safe DNA Gel Stain
- Promega Wizard® SV Gel and PCR Clean-Up System

## Procedure

1- Add to a sterile 1.5mL tube on ice:

Component	72 $\mu$ L reaction (GFP)	60 $\mu$ L reaction (variable part of the toehold switch)
Insert	60 $\mu$ L	50 $\mu$ L
Buffer	rCutSmart = 7.2 $\mu$ L	3.1 buffer = 6 $\mu$ L
Restriction enzymes	XhoI = 0.6 $\mu$ L EcoRI = 1.2 $\mu$ L	BglII = 1 $\mu$ L BamHI = 1.5 $\mu$ L
Nuclease-free H <sub>2</sub> O	3 $\mu$ L	1.5 $\mu$ L

2- Incubate 12 hours at 37°C.

3- Purify your digested insert via gel migration. Prepare a 1.25% low melting gel, run your sample for 45 minutes et 25V.

4- Under UV, cut the insert band. Weight your gel slice using a precision scale.



5- Perform a PCR clean-up using the Promega Wizard® SV Gel and PCR Clean-Up System.

Add  $x \mu\text{L}$  of Membrane Binding Solution per  $x \text{ mg}$  of gel to your gel slices. Heat your sample at  $50^{\circ}\text{C}$  for at least 10 minutes and follow the 5.A. DNA Purification by Centrifugation part of the protocol given in the Promega Wizard® SV Gel and PCR Clean-Up System kit.



# Plasmid digestion

## Introduction

This is the protocol for plasmid pET14d(+) digestion using restriction enzymes for molecular cloning.

## Material

- Ice bucket
- DNA insert (~100-200 ng/ $\mu$ L)
- NEB® rCutSmart Buffer
- NEB® 3.1 Buffer
- NEB® XhoI
- NEB® XhoI
- NEB® XhoI
- NEB® EcoRI
- Nuclease-free water
- TAE 1X
- Low Melting Agarose powder
- SYBR™ Safe DNA Gel Stain
- Promega Wizard® SV Gel and PCR Clean-Up System

## Procedure

1- Add to a sterile 1.5mL tube on ice:

component	60 $\mu$ L reaction (cloning with GFP)	60 $\mu$ L reaction (cloning with variable part of the toehold switch)
DNA plasmid	50 $\mu$ L	50 $\mu$ L
Buffer	rCutSmart = 6 $\mu$ L	3.1 Buffer = 6 $\mu$ L
Restriction enzymes	EcoRI = 1.5 $\mu$ L XhoI = 1 $\mu$ L	BglII = 1 $\mu$ L BamHI = 1.5 $\mu$ L
Nuclease-free H <sub>2</sub> O	1.5 $\mu$ L	1.5 $\mu$ L

2- Incubate 12 hours at 37°C.

4- Under UV, cut the upper band (open plasmid). Weight your gel slice using a precision scale.

5- Perform a PCR clean-up using the Promega Wizard® SV Gel and PCR Clean-Up System.

Add x  $\mu$ L of Membrane Binding Solution per x mg of gel to your gel slices. Heat your sample at 50°C for at least 10 minutes and follow the 5.A. DNA Purification by



Centrifugation part of the protocol given in the Promega Wizard® SV Gel and PCR Clean-Up System kit.

## Ligation

### Introduction

This is the protocol for the ligation of the digested plasmid with the digested insert using NEB® T4 DNA ligase.

### Material

- NEB® T4 DNA ligase
- NEB® 1X T4 DNA Ligase Reaction Buffer
- Digested plasmid (20ng/μL)
- Digested insert (40ng/μL)

### Procedure

1- Add to a sterile 1.5mL tube:

Component	20 μL reaction
Digested plasmid	4 μL
Digested insert	5 μL
NEB® 1X T4 DNA Ligase Reaction Buffer	10 μL
NEB® T4 DNA ligase	1 μL

2- Leave at room temperature for 10 minutes

3- Store at -20°C



# BI21(DE3) E. coli heat shock transformation

## Introduction

This is the protocol for the heat shock transformation of BI21(DE3) E. coli with your plasmid assembled using the golden gate method previously detailed.

## Material

- Bucket of ice
- BI21(DE3) E. coli (120  $\mu$ L)
- Plasmid DNA (~100 ng/ $\mu$ L)
- Terrific Broth (TB) media
- Round-bottom tube
- Petri dishes (LB + antibiotic of your plasmid resistance cassette)

## Procedure

1- Store all your material on ice under the hood. Bacteria should thaw slowly (~10 minutes).

At the same time, prepare aliquots of TB and incubate them at 37°C. Turn on your water bath at 42°C.

2- Add 2  $\mu$ L of the plasmid at the bottom of the P20 pipette

3- Add 120  $\mu$ L of the bacteria, gently tap the tube to mix.

4- Leave to incubate on ice for 30 minutes. Gently tap the tube every 10 minutes.

5- Heat shock 30 seconds in a water bath at 42°C

6- Leave 2 minutes on ice

7- Add 150  $\mu$ L of the TB media

8- Incubate on the shaker 220 RPM, 37°C for 90 minutes

9- Plate 120  $\mu$ L on your petri dish

10- Incubate at 37°C overnight

