

3. PCR & Agarose Gel Electrophoresis

Material

2 x Phanta Max Master Mix
Primer
Double Distillation Water (ddH₂O)
Loading Buffer
DNA Marker

Step

- ① Adding 25 μ l of 2 x Phanta Max Master Mix, <30ng Template DNA, 4 μ l of primer mix, then add ddH₂O up to 50 μ l into a PCR tube. Different group use specific tube with distinctive sign.
- ② Place those PCR tubes into Peltier thermal cycler.
- ③ Set the protocol as follow: begin at 95 $^{\circ}$ C for 30 secs, then keep 95 $^{\circ}$ C for 15 secs for denaturation, decrease to 60 $^{\circ}$ C for 15 secs, 72 $^{\circ}$ C for 1 min/kb and repeat that cycle 30 times, finally maintain 16 $^{\circ}$ C infinity.
- ④ Adding all the samples to the hole. Run the gel at 120V for about 25 minutes. Check the result under the Blue Light Gel Imager.
- ⑤ *Note*
- ⑥ As we use Green Taq Mix, we needn't add Loading when Agarose Gel Electrophoresis.