

16. Protocol for Western blotting

- ① Preparation of lysate from 60 mm×15 mm cell culture plate:
 - a. Aspirate the growth medium from the cells, and wash the cells with PBS.
 - b. Aspirate the PBS, then add 600 μ L trypsin EDTA solution.
 - c. Incubation for several minutes at 37°C and check the culture with inverted microscope to be sure that the cells are rounded up which indicates they are detached from the surface.
 - d. Add 2 mL fresh complete medium, then pipet all the adherent cells into cell suspension and pipet up and down.
 - e. Transfer the cell suspension into a pre-cooled 15 mL centrifuge tube and centrifuge at 1200 rpm for 5 minutes.
 - f. Aspirate the filtrate, add 1mL PBS and resuspend the cell, transfer the cell suspension into a pre-cooled 1.5 mL microcentrifuge tube.
 - g. Centrifuge at 1200 rpm for 5 minutes at 4°C, then aspirate the filtrate.
 - h. Add 35 μ L RIPA Lysis Buffer containing protease inhibitor cocktail and phospho-stop, and mix or vortex briefly to resuspend the cells completely
 - i. Incubate on ice for 15 minutes. Vortex briefly to resuspend and lyse residual cells.
 - j. Centrifuge at 13200 rpm for 15 minutes at 4°C.
 - k. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice, and discard the pellet.
- ② Protein Quantification:

We use Beyotime BCA Protein Assay Kit.

 - a. Dilute 10 μ L BSA Standard (5 mg/mL) to 100 μ L.
 - b. Prepare Working Reagent by mixing 50 parts of BCA™ Reagent A with 1 part of BCA™ Reagent B (50:1, Reagent A:B).
 - c. i. Pipette 0、1、2、4、6、8、12、16、20 μ L of standard into a 96 microplate well and dilute each standard to 20 μ L. ii. Pipette 20 μ L of each diluted sample replicate into the microplate well. iii. Add 200 μ L of the Working Reagent to each well and mix plate thoroughly by pipette. iv. Cover plate and incubate at 37°C for 30 minutes.
 - d. Measure the absorbance at 595 nm on a plate reader.
 - e. Prepare a standard curve by plotting the average Blank-corrected 595 nm measurement for each BSA standard vs. its concentration in μ g/ml. Use the standard curve to determine the protein concentration of each unknown sample.
- ③ Reduce and denature samples:

Boil each cell lysate in LDS sample buffer(4 \times) at 95°C for 5 minutes. Lysates can be aliquoted and stored at -20°C for future use.
- ④ Loading and running the gel:

Load 60–80 μ g of total protein from cell lysate into the wells of the SDS-PAGE gel(10% Gel percentage), along with pre-stained molecular weight marker.
- ⑤ Run the gel for 55 minutes at 180 V.
- ⑥ Transfer proteins from the gel to membrane:
 - a. i. Prepare transfer buffer(1L): 3.03 g Trizma base (25 mM), 14.4 g Glycine (192 mM), 200 mL methanol(20%), 800 mL water; ii. Prepare PBST(1L):PBS(1L),

0.5ml Tween 20 (0.05%)

- b. Cut a piece of PVDF membrane according to the size of gel and wet in methanol for 2 minutes.
- c. Transfer the membrane to the transfer buffer.
- d. Assemble transfer stack as follows:



- e. Attach the electrodes. Set the power supply to 350 mA (constant current) for 1h 20 min at 4°C.

⑦ Blocking:

- a. Prepare Blocking buffer: Add 2.5g non-fat powdered milk (5%) to 50 mL PBST buffer, mix thoroughly.
- b. Remove the membrane to a small container and stain in Ponceau stain reagent for several minutes, then wash for short time with PBST to check the transfer.
- c. Cut the membrane into several stripes according to the size of target protein.
- d. Completely destain in PBST.
- e. Incubate the membrane in 5% milk for 1 h at room temperature.

⑧ Primary antibody incubation:

Incubate the membrane with appropriate dilutions of primary antibody in blocking buffer, 4°C overnight at recommended concentration

⑨ Secondary antibody incubation:

- a. Wash the membrane in 5 washes of PBST, 6 min each.
- b. Incubate the membrane with the recommended dilution of conjugated secondary antibody in blocking buffer at room temperature for 1 h.

⑩ ECL detection: We use Pierce™ ECL Western Blotting Substrate.

- a. Wash the membrane in 3 washes of PBST, 5 min each.
- b. Prepare the substrate working solution by mixing equal parts of Detection Reagents 1 and 2. Use 0.125mL Working Solution per cm² of membrane.
- c. Incubate blot with working solution for 1 minute at RT.
- d. Remove blot from working solution and place it in a clear plastic wrap. Use an absorbent tissue to remove excess liquid.
- e. Place the membrane in a chemiluminescence imaging system with the protein side facing up and carefully press out any bubbles from between the chemiluminescence imaging system and the membrane.
- f. Cover the chemiluminescence imaging system and detect