

03_04 PCR for ssDNA Samples

MITTWOCH, 2.6.2021

Goal setting

- Performing a PCR reaction using the Phusion™ High-Fidelity DNA Polymerase

Terms / abbreviations

- DMSO = Dimethyl sulfoxide
- PCR = Polymerase chain reaction

Risk areas

- If spilled, always wipe surface with alcohol

 Hazard symbols



Required materials/ information

- Chemicals:
 - 5x Phusion HF buffer, ThermoFisher
 - 5x Phusion GC buffer, ThermoFisher
 - 10 mM dNTP, ThermoFisher
 - 10 µM Forward primer, Ella Biotech
 - 10 µM Reverse primer, Ella Biotech
 - 50 mM MgCl₂ solution, ThermoFisher
 - Autoclaved MilliQ water, Sartorius arium pro VF
 - DMSO, ThermoFisher
 - Phusion DNA Polymerase (2 U/µL), ThermoFisher
 - Template DNA
- Material:
 - Ice or cooling rack
 - PCR tubes, Sarstedt
 - Trash bags, Th. Geyer GmbH & Co. KG

Templates, devices, software

- Nanodrop spectrophotometer, ThermoFisher
- Thermocycler, Eppendorf
- Pipettes, Eppendorf

Preliminary work

- Any reaction creating DNA that should be amplified
- [04_01 Spectralphotometer NanoDrop](#) to detected desired amount of template DNA

Operation

- Put the PCR tube on ice or the cooling rack
- Add the components in the order listed in the following table

Composition of Master Mix					
	A	B	C	D	E
1	Component	20 [μL]	50 [μL]	Final concentration	Unit
2	MilliQ water	add to 20 μL	add to 50 μL		
3	5X Phusion HF Buffer	4	10	1	x
4	10 mM dNTPs	0.4	1	200	μM
5	Forward primer 10 μM	1	2.5	0.5	μM
6	Reverse primer 10 μM	1	2.5	0.5	μM
7	Template DNA *	x	x	0.2	
8	Phusion High-Fidelity DNA Polymerase	0.2	0.5	0.02	U/μl

* use 1 μL from standard TdT reaction as template

- Run the PCR program in the thermocycler according to the following table

PCR Program for Thermocycler				
	A	B	C	D
1	Temperature [°C]	Time	Unit	Repeats
2	98	30	s	1x
3	98	5	s	30x
4	63	10	s	
5	72	30	s	
6	72	5	min	1x

Increase elongation time for fragments > 1 kb

Disposal

- Autoclave trash bags, discard in S1 waste

Troubleshooting

- Wear gloves to reduce the risk of DNase and RNase contamination
- The Phusion DNA Polymerase should be pipetted carefully and gently as the high glycerol content (50%) in the storage buffer may otherwise lead to pipetting errors

Sources

 MAN0012393_Phusion_HighFidelity_DNAPolymerase_UG.pdf

Follow-up work

- Purify fragments with a cleanup-kit
- If materials are empty, care about new order