

Plasmid Preparation Procedure (By Hand)

Introduction

Plasmid Preparation Procedure (by hand) is used only when necessary.

Materials

› Solutions

- › TENS (make fresh each time) To 4.5 mL TE add: 250 λ 10% SDS 250 λ 2N NaOH
- › Sodium Acetate 3.0M; pH 5.5
- › TE with 10ug/mL of RNase A

Procedure

Main Procedure

- ✓ 1. Spin 1.5 mL of overnight culture for 30 seconds in microfuge.
- ✓ 2. Aspirate off all but 100 λ of the supernatant and resuspend the pellet by vortexing.
- ✓ 3. Add 300 λ of TENS and mix by inversion for 3-5 minutes. The solution should become viscous.
- ✓ 4. Add 150 λ of sodium acetate and vortex. A fine white precipitate should form.
- ✓ 5. Centrifuge for 2.5 minutes at 10K.
- ✓ 6. TRANSFER the supernatant to a clean tube and add 2 volumes (1 mL) of room temperature EtOH.
- ✓ 7. Vortex and pellet DNA by centrifugation for 2-5 minutes at 10K.
- ✓ 8. Wash pellet with 70% ethanol and allow the pellet to dry.
- ✓ 9. Resuspend the pellet in 30 λ of TE with RNaseA (with a P100).
- ✓ 10. Digest 5-10 λ as usual.