



# AMALTHEA



## PROTOCOLS

COMMONLY USED  
REAGENTS RECIPES

## Commonly Used Reagents Recipes

### Acid precipitation solution

- 1 M HCl
- 0.1 M sodium pyrophosphate

*Nucleic acids can also be precipitated with a 10% (w/v) solution of trichloroacetic acid (TCA); however, this recipe is cheaper, easier to prepare, and just as efficient.*

### Ammonium acetate, 10 M

- Dissolve 385.4 g ammonium acetate in 150 ml H<sub>2</sub>O
- Add H<sub>2</sub>O to 500 ml

### BBS (BES-buffered solution), 2×

- 50 mM N,N-bis(2-hydroxyethyl)-2-aminoethanesulfonic acid (BES; Calbiochem)
- 280 mM NaCl
- 1.5 mM Na<sub>2</sub>HPO<sub>4</sub>, pH 6.95
- 800 ml H<sub>2</sub>O
- Adjust pH to 6.95 with room temperature 1 M NaOH
- H<sub>2</sub>O to 1 liter

Filter sterilize through a 0.45-μm nitrocellulose filter (Nalgene)

Store in aliquots at −20°C (can be frozen and thawed repeatedly)

*The pH of this solution is critical (pH 6.95 to 6.98). When a new batch of 2× BES buffer is prepared, its pH should be checked against a reference stock prepared (and tested) earlier.*

### CaCl<sub>2</sub>, 1 M

- 147 g CaCl<sub>2</sub>·2H<sub>2</sub>O
- H<sub>2</sub>O to 1 liter

**Denhardt solution, 100×**

- 10 g Ficoll 400
- 10 g polyvinylpyrrolidone
- 10 g bovine serum albumin (Pentax Fraction V; Miles Laboratories)
- H<sub>2</sub>O to 500 ml

*Filter sterilize and store at -20°C in 25-ml aliquots*

**Dithiothreitol (DTT), 1 M**

- Dissolve 15.45 g DTT in 100 ml H<sub>2</sub>O
- Store at -20°C

**EDTA (ethylenediamine tetraacetic acid), 0.5 M (pH 8.0)**

- Dissolve 186.1 g Na<sub>2</sub>EDTA·2H<sub>2</sub>O in 700 ml H<sub>2</sub>O
- Adjust pH to 8.0 with 10 M NaOH (~50 ml)
- Add H<sub>2</sub>O to 1 liter

**HBSS (Hanks balanced salt solution)**

- 5.4 mM KCl
- 0.3 mM Na<sub>2</sub>HPO<sub>4</sub>
- 0.4 mM KH<sub>2</sub>PO<sub>4</sub>
- 4.2 mM NaHCO<sub>3</sub>
- 1.3 mM CaCl<sub>2</sub>
- 0.5 mM MgCl<sub>2</sub>
- 0.6 mM MgSO<sub>4</sub>
- 137 mM NaCl
- 5.6 mM D-glucose
- 0.02% phenol red (optional)
- Add H<sub>2</sub>O to 1 liter and adjust pH to 7.4

*HBSS can be purchased from Biofluids or Whittaker.*

*HBSS may be made or purchased without CaCl<sub>2</sub> and MgCl<sub>2</sub>. These are optional components that usually have no effect on an experiment. In some cases, however, their presence may be detrimental to a procedure. Consult the individual protocol to see if the presence or absence of these components is recommended in the materials list.*

### **HCl, 1 M**

Mix in the following order:

- 913.8 ml H<sub>2</sub>O
- 86.2 ml concentrated HCl

### **HeBS (HEPES-buffered saline) solution, 2×**

- 16.4 g NaCl
- 11.9 g HEPES acid
- 0.21 g Na<sub>2</sub>HPO<sub>4</sub>
- 800 ml H<sub>2</sub>O
- Titrate to pH 7.05 with 5 M NaOH
- Add H<sub>2</sub>O to 1 liter

Filter sterilize through a 0.45-μm nitrocellulose filter

Test for transfection efficiency and store at –20°C in 50-ml aliquots

*An exact pH is extremely important for efficient transfection. The optimal pH range is 7.05 to 7.12.*

### **KCl, 1 M**

- 74.6 g KCl
- H<sub>2</sub>O to 1 liter

### **MgCl<sub>2</sub>, 1 M**

- 20.3 g MgCl<sub>2</sub>·6H<sub>2</sub>O
- H<sub>2</sub>O to 100 ml

### **MgSO<sub>4</sub>, 1 M**

- 24.6 g MgSO<sub>4</sub>·7H<sub>2</sub>O
- H<sub>2</sub>O to 100 ml

**MOPS buffer**

- 0.2 M MOPS [3-(N-morpholino)propanesulfonic acid], pH 7.0
- 0.5 M sodium acetate
- 0.1 M EDTA

Store in the dark and discard if it turns yellow.

**NaCl, 5 M**

- 292 g NaCl
- H<sub>2</sub>O to 1 liter

**NaOH, 10 M**

- Dissolve 400 g NaOH in 450 ml H<sub>2</sub>O
- Add H<sub>2</sub>O to 1 liter

**PBS (phosphate-buffered saline)**

10× stock solution, 1 liter:

- 80 g NaCl
- 2 g KCl
- 11.5 g Na<sub>2</sub>HPO<sub>4</sub>·7H<sub>2</sub>O
- 2 g KH<sub>2</sub>PO<sub>4</sub>

Working solution, pH ~7.3:

- 137 mM NaCl
- 2.7 mM KCl
- 4.3 mM Na<sub>2</sub>HPO<sub>4</sub>·7H<sub>2</sub>O
- 1.4 mM KH<sub>2</sub>PO<sub>4</sub>

**Potassium acetate buffer, 0.1 M**

- Solution A: 11.55 ml glacial acetic acid/liter (0.2 M).
- Solution B: 19.6 g potassium acetate (KC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>)/liter (0.2 M).

Referring to Table A.2.2 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H<sub>2</sub>O to 100 ml. *This may be made as a 5- or 10-fold concentrate by scaling up the amount of potassium acetate in the same volume. Acetate buffers show concentration-dependent pH changes, so check concentrate pH by diluting an aliquot to the final concentration. To prepare buffers with pH intermediate between the points listed in Table A.2.2, prepare closest higher pH, then titrate with solution A.*

**Table A.2.2** Preparation of 0.1 M Sodium and Potassium Acetate Buffers<sup>a</sup>

Desired pH	Solution A (ml)	Solution B (ml)
3.6	46.3	3.7
3.8	44.0	6.0
4.0	41.0	9.0
4.2	36.8	13.2
4.4	30.5	19.5
4.6	25.5	24.5
4.8	20.0	30.0
5.0	14.8	35.2
5.2	10.5	39.5
5.4	8.8	41.2
5.6	4.8	45.2

<sup>a</sup>Adapted by permission from CRC (1975).

#### Potassium phosphate buffer, 0.1 M

- Solution A: 27.2 g KH<sub>2</sub>PO<sub>4</sub> per liter (0.2 M).
- Solution B: 34.8 g K<sub>2</sub>HPO<sub>4</sub> per liter (0.2 M).

Referring to Table A.2.3 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H<sub>2</sub>O to 200 ml. *This may be made as a 5- or 10-fold concentrate by scaling up the amount of potassium phosphate in the same volume. Phosphate buffers show concentration-dependent pH changes, so check concentrate pH by diluting an aliquot to the final concentration.*

**Table A.2.3** Preparation of 0.1 M Sodium and Potassium Phosphate Buffers<sup>a</sup>

Desired pH	Solution A (ml)	Solution B (ml)	Desired pH	Solution A (ml)	Solution B (ml)
5.7	93.5	6.5	6.9	45.0	55.0
5.8	92.0	8.0	7.0	39.0	61.0
5.9	90.0	10.0	7.1	33.0	67.0
6.0	87.7	12.3	7.2	28.0	72.0
6.1	85.0	15.0	7.3	23.0	77.0
6.2	81.5	18.5	7.4	19.0	81.0
6.3	77.5	22.5	7.5	16.0	84.0
6.4	73.5	26.5	7.6	13.0	87.0
6.5	68.5	31.5	7.7	10.5	90.5
6.6	62.5	37.5	7.8	8.5	91.5
6.7	56.5	43.5	7.9	7.0	93.0
6.8	51.0	49.0	8.0	5.3	94.7

<sup>a</sup>Adapted by permission from CRC (1975).

**SDS electrophoresis buffer, 5×**

- 15.1 g Tris base
- 72.0 g glycine
- 5.0 g SDS
- H<sub>2</sub>O to 1000 ml
- Dilute to 1× or 2× for working solution, as appropriate

Store up to 1 month at 0° to 4°C

*Do not adjust the pH of the stock solution, as the solution is pH 8.3 when diluted.*

**SED (standard enzyme diluent)**

- 20 mM Tris·Cl, pH 7.5
- 500 µg/ml bovine serum albumin (Pentax Fraction V)
- 10 mM 2-mercaptoethanol

Store up to 1 month at 4°C

**Sodium acetate, 3 M**

- Dissolve 408 g sodium acetate·3H<sub>2</sub>O in 800 ml H<sub>2</sub>O
- Add H<sub>2</sub>O to 1 liter

Adjust pH to 4.8 or 5.2 (as desired) with 3 M acetic acid

**Sodium acetate buffer, 0.1 M**

**Solution A: 11.55 ml glacial acetic acid/liter (0.2 M).**

**Solution B: 27.2 g sodium acetate (NaC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>·3H<sub>2</sub>O)/liter (0.2 M).**

Referring to Table A.2.2 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H<sub>2</sub>O to 100 ml. (See Potassium acetate buffer recipe for further details.)

**Sodium phosphate buffer, 0.1 M**

- Solution A: 27.6 g NaH<sub>2</sub>PO<sub>4</sub>·H<sub>2</sub>O per liter (0.2 M).
- Solution B: 53.65 g Na<sub>2</sub>HPO<sub>4</sub>·7H<sub>2</sub>O per liter (0.2 M).

Referring to Table A.2.3 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H<sub>2</sub>O to 200 ml. (See Potassium phosphate buffer recipe for further details.)

**SSC (sodium chloride/sodium citrate), 20×**

- 3 M NaCl (175 g/liter)
- 0.3 M Na<sub>3</sub>citrate·2H<sub>2</sub>O (88 g/liter)
- Adjust pH to 7.0 with 1 M HCl

**STE buffer**

- 10 mM Tris·Cl, pH 7.5
- 10 mM NaCl
- 1 mM EDTA, pH 8.0

**TE (Tris/EDTA) buffer**

- 10 mM Tris·Cl, pH 7.4, 7.5, or 8.0 (or other pH; see recipe)
- 1 mM EDTA, pH 8.0

**TEA (triethanolamine) solution**

- 50 mM triethanolamine, pH ~11.5
- 0.1% Triton X-100
- 0.15 M NaCl

*Add Triton X-100 as a 10% stock sterilized by Millipore filtration and stored in the dark to prevent photooxidation (stock is stable 5 years at room temperature).*

**TEN (Tris/EDTA/NaCl) solution**

- 40 mM Tris·Cl, pH 7.5
- 1 mM EDTA, pH 8.0
- 150 mM NaCl

Store up to 6 months at room temperature

**TM buffer, 10×**

- 100 mM Tris·Cl, pH 8.0
- 100 mM MgCl<sub>2</sub>



### **Tris-buffered saline (TBS)**

- 100 mM Tris·Cl, pH 7.5
- 0.9% (150 mM) NaCl

Store up to several months at 4°C

### **Tris·Cl [tris(hydroxymethyl)aminomethane], 1 M**

- Dissolve 121 g Tris base in 800 ml H<sub>2</sub>O
- Adjust to desired pH with concentrated HCl
- Mix and add H<sub>2</sub>O to 1 liter

*Approximately 70 ml of HCl is needed to achieve a pH 7.4 solution, and approximately 42 ml for a solution that is pH 8.0.*

*IMPORTANT NOTE: The pH of Tris buffers changes significantly with temperature, decreasing approximately 0.028 pH units per 1°C. Tris-buffered solutions should be adjusted to the desired pH at the temperature at which they will be used. Because the pK<sub>a</sub> of Tris is 8.08, Tris should not be used as a buffer below pH ~7.2 or above pH ~9.0.*

### **TTBS (Tween 20/TBS)**

- 0.1% Tween 20 in Tris-buffered saline (TBS; see recipe)

Store up to several months at 4°C

### **1 M IPTG (1 M IPTG Recipe) Stock**

1. Dissolve 2.38 g of IPTG in 8 mL of distilled H<sub>2</sub>O.
2. Bring to a final volume of 10 mL with molecular biology grade H<sub>2</sub>O.
3. Filter sterilize with a 0.22 µ syringe filter.
4. Store in 1mL aliquots at -20 °C.



# AMALTHEA

## PROTOCOLS

