

Colony PCR

MARDI 17/08/2021

Preparing primers:

180 μ L RNA free water + 20 μ L of primer stock (->10 fold dilution)

Setting up the reactions:

- Collect a single colony into 50 μ L of RNase free water and lyse at 95°C for 2-5 minutes.
- Set up the reactions like so:

Table1			
	A	B	
1	Reagent	Volume	
2		1x	
3	Nuclease-free water	19 μ L	
4	2x DreamTaq Master Mix	25 μ L	
5	Forward Primer (10 μ M)	2.5 μ L	
6	Reverse Primer (10 μ M)	2.5 μ L	
7	Template	1 μ L	
8	Total Volume	50 μ L	

*A master mix with the all components except for the template can be made to limit pipetting.
Vortex on very low to not damage the enzymes*

- Set up the thermal cycler like so:

Table2



	A	B	C	D	E
1	Thermocycler Protocol: DreamTaq Master Mix				
2		Temp	Time		
3	Start	95 °C	1 min	Melt	
4					
5	Cycle 1	95 °C	30 sec	Melt	35 Cycles
6	Cycle 2	depends on annealing temps of primers*	30 sec	Anneal	
7	Cycle 3	72 °C	1 min	Extend	
8					
9	Finish	72 °C	5 min	Extend	
10	Store	10 °C	Forever	Store	

*55°C is a standard temperature