

Washing and preparing *K. Phaffii* cryo stocks ready for transformation

Introduction

As our previous experiments have proved to show, *K. phaffii* cells can be washed and prepared for electroporation, and be put in the freezer for storage until the electroporation can occur. This is a great way of saving time, since the washing step itself takes at least half a day.

This protocol is based on a protocol adapted from Adrian and provided by Louise, called "Transformation into *K. phaffii*: Electroporation"

Materials

› Pipettes

- › P5000
- › P1000

› Reagents

- › YPD
- › 100 mM LiAc, 10 mM DTT, 0.6 M sorbitol, and 6 mM Tris HCl pH 8.5
- › 1M sorbitol

Procedure

Preparation of *K. phaffii* cultures

1. Inoculate 4 mL of YPD with a single cell colony of *K. phaffii* GS115 Δ KU70 and incubate O/N at 30°C

Note: The inoculation should be done early afternoon, since the Δ KU70 strains grows slower than a wild type

2. Transfer 500 μ L from the pre-culture to a 50mL falcon tubes with 25mL YPD
3. Incubate for 5 hours at 30°C and 200RPM

Preparation of Master Mix for the washing

4. Mix the necessary amount according to the table below

Master Mix recipe							
	A	B	C	D	E	F	G
1	Reagent	Volume for 500 mL Mastermix [mL]	Volume per transformation reaction [μL]	Mastermix volume for 8x8mL [μL]			
2	Sorbitol 1M	300	4800	42240		Number of reactions	8
3	LiAc 1 M	50	800	7040		Safety margin (%)	10
4	Tris HCl pH 8.5 1M	3	48	422.4			
5	DTT	5	80	704			
6	sterile dH2O	142	2272	19993.6			
7	Total volume [mL]	500	8 mL	70.4 mL			

5. CRITICAL ADD THE DTT JUST BEFORE ADDING THE MASTER MIX TO THE CELL PELLETS

Washing of *K. phaffii* cells

6. Pellet cells at 3500G for 5 minutes, and discard the supernatant

7. Gently resuspend cells in 8mL Master Mix (of 100 mM LiAc, 10 mM DTT, 0.6 M sorbitol, 6 mM Tris HCl pH 8.5), then transfer to new 50mL falcon tubes, and incubate at room temperature for 30 minutes.

Make sure cells do not precipitate, gently swirl every 10 minutes.

8. Pellet cells at 3500G for 5 minutes, and discard supernatant.

9. **Work on ice:** Resuspended cells in 1.5mL 1M ice-cold sorbitol and transfer to sterile 2mL eppendorff tubes.

10. **Work on ice:** Wash cells 3 times by centrifugation at 13000G for 30 seconds by resuspending in 1ml 1M ice-cold sorbitol, vortex between washes. After 3rd time add 1ml 1M ice-cold sorbitol.

11. Work on ice: Distribute cells in 5 sterile 1.5mL cryo tubes (200μL in each)

12. Immediately put the cryo tubes in the -80°C freezer

Keep on ice until putting them in the freezer

Note: It is a good idea to also have pre-cooled the cryo tubes before adding the washed cells

1 tube = 1 electroporation reaction