

12. Protocol for SEAP assay in vitro

Material

Chemiluminescence-based SEAP assay kit (e.g., Sigma Aldrich; cat. no. 1177984200)
buffer (2x): 20 mM L-Homoarginine hydrochloride, 1 mM MgCl₂, 21% (v/v)
diethanolamine, pH 9.8.

Substrate for SEAP: 120 mM para-nitrophenyl phosphate (pNPP) in 2SEAP buffer.

Transparent 96-well plates.

Colorimetric plate-reader.

Steps

- ① Sample 200 μ L culture medium from each well, heat inactivate at 65 C for 30 min.
- ② During the heat inactivation procedure, warm up 2 SEAP buffer (100 μ L/well) at 37 °C.
- ③ Add 1/5 buffer volume of pNPP (20 μ L/well) substrate into the 2x buffer to prepare the "Detection Mixture."
- ④ Add 80 μ L heated medium into the 96-well plate, add 120 μ L Detection Mixture.
- ⑤ Measure absorption at 405 nm, 30 s per read for 10 reads.
- ⑥ Calculate enzymatic activity.

Note

SEAP levels in culture supernatants can be quantified with a pnitrophenylphosphate-based colorimetric assay