

## Transfection

1. 100  $\mu$ l of the cell culture ( **$10^8$**  cells per ml) is fixed with 200  $\mu$ l fixing solution
2. 10  $\mu$ l of the cell culture are placed in a counting chamber
3. The cells are counted under a microscope
4. The cells are centrifuged at room temperature (RT) by **3310 rpm for 5 min**
5. Take the supernatant off and wash the pellet with **1 ml of transfection buffer**
6. The cells are centrifuged again at RT by **3310 rpm for 5 min** and take the supernatant off
7. Add 150  $\mu$ l transfection buffer per transfection to the pellet
8. Before the DNA constructs are added to the cell culture, they are sterilized for **5 min by 95°C** in the PCR cycler
9. Put **150  $\mu$ l cell culture** in electroporation cuvettes and add **1–20  $\mu$ g of the DNA (~50  $\mu$ l)** construct.
10. For a **negative control**, add 50  $\mu$ l of distilled water to the 150  $\mu$ l cell culture in the electroporation cuvettes
11. The cuvette with the cell / DNA suspension is placed in the Nucleofector and the transfection is carried out with the **program U-033**
12. The suspension is pipetted from the cuvette into **1 ml BHI medium** ( → **Basic Protocol**: BHI medium for the cellcultures aftee transfection) into an Eppi and **incubated at 27 ° C** overnight