

Cell Culture/Cell Passage

Introduction

This protocol will help you to know which material you need and how do you process to successfully culture your cell

Material

- 500ml of DMEM (Invitrogen)
- 50 ml of FBS
- HeLa cells or MCF-7 cells
- T175 flasks
- Gibco Trypsin-EDTA (0,5%)
- PBS

Culture Cell :

- 1- Prepare media by Mixing the 500ml of DMSO and 50mL of FBS and filtering it.
- 2- Count your cells to plan how many cells you will initially put inside the flask.
- 3- add the desired volume in 2ml of media. The number of cells needed may differ with the cell line you are working with.
- 4- Put 28 mL of media inside the T175 flask
- 5- Put your cell + media inside the T175 Flask. And put the flask in the incubator at 37°C and 5% of CO₂

Cell Passage

Prerequisite : Cell at 80% of confluency

- 1- Aspirate the medium (put your Pasteur tip on one of the top corners of the flask so you don't touch the cells)



- 2- Add 5ml of PBS to wash the cells and remove the dead ones (move laterally the flask to wash all the cells)
- 3- Aspirate the PBS
- 4- Add 5ml of PBs
- 5- Aspirate
- 6- Add 3mL of PBS and 1mL of trypsin
- 7- Incubate at 37°C for 5-10min
- 8- Go back under the hood: tap the flask to be sure you've detached all the cells
- 9- Add 4mL of medium, with your 10mL pipette, wash the flask and homogenise the solution (DMEM stops the trypsin action). Avoid bubbles as much as you can
- 10- Transfer the 8mL in a 15mL Falcon tube
- 11- Spin for 4min at 1200rpm, 20°C
- 12- Aspirate the supernatant
- 13- Tap the pellet to break it and disperse the cells
- 14- Add 1mL medium with 5mL pipette
- 15- Resuspend by adding 9ml of medium
- 16- Count the cells:
- 17- Label your T175 flask
- 18- Add 29.5 ml of medium + 550 μ L of cells

