


USER Cloning


Goal: To scarlessly combine two linear fragments of DNA, each of which contain uracil residues, into a larger linear or circular fragment.

See figures below that explain the concept behind USER cloning (for a hypothetical situation where one wants to scarlessly ligate Part1 to Part2)


PCR **Part 1** off of Template 1



PCR Part 2 off of Template 2



Digest PCR Products with DpnI and User Enzyme



Protocol:

(1) Amplify target sequences (Parts 1 and 2) with primers that contain an adenine at the 5' end and a uracil roughly 8-11 bp later (in place of a thymine). In between the 5' adenine and the uracil will be the USER cloning overhang sequence - which may be standardized or part of the gene you are trying to clone. When designing overhangs, I try to find an A and T roughly 8-11 bp away from each other with a lower GC content (higher specificity). Use Pfu Turbo Cx polymerase to minimize PCR errors - although I have only had success with this polymerase for amplicons \leq 2 kb.

Pfu Turbo Cx Amplification Reaction Mix (20 μ L)

1 μ L Pfu Turbo Cx
1 μ L 10 μ M Fwd Primer
1 μ L 10 μ M Rev Primer
1 μ L DNA template
1 μ L 10 mM dNTPs
2 μ L Pfu Turbo Cx Buffer (10x)
13 μ L H₂O

Thermocycler Settings:

95°C for 3 minutes (or 5 minutes for colony/genomic DNA templates)

30 Cycles of:

95°C for 30 seconds

55°C for 30 seconds

68°C for extension time (1 kb/minute) + some (I'm generous here)

1 Cycle of:

68°C for 20 minutes

Hold at 4°C

(2) Perform USER cloning reaction:

Note, run reactions with just each part on its own as negative controls to see how much "background" you are getting from each one - replace missing parts with H₂O. This is most important for parts that contain an origin of replication and an antibiotic resistance marker

2-way assembly Reaction Mixture (4.8 uL)

0.4 uL DpnI
0.4 uL USER enzyme
2 uL Part 1
2 uL Part 2

3-way assembly Reaction Mixture (4.8 uL)

0.4 uL DpnI
0.4 uL USER enzyme
1.33 uL Part 1
1.33 uL Part 2
1.33 uL Part 3

On the thermocycler:

37°C for 1 hour.
22°C (room temperature) for 15 minutes.
Cool to 4°C.

Directly transform 1 uL of the product into chemically competent cells (follow [plasmid transformation](#) protocol) - Plate everything and should get 20 - 50 colonies.

Credit: Spencer Glantz