



Look What They've Done To My Shoes!

SCU_China Project 2016

Protocol
AarC Team



AarC

Overlap extension PCR Cecropin XJ Expression Team

1. Product two fragment

Primers

(1) Previous fragment:

gctctagatagctactagagtacacaggaagtagtagtgaccgaacgtatc Xba-I-AarC-F
gcacgatcgaaatagctcttgacgatcggacggtaatccgg Repair-Pst-I-R

(2) Rearward fragment:

ccggattaccgtccgatgctgcaagactatttcgatcgtgc Repair-Pst-I-F
tgactgcagcggcccgtactagtaaaaaaccctcaagaccgtagagg T7T-Pst-I-R

Diluted to 10 pmol/ul

Reaction system

(1) 5ul of I-5™ 2X High-Fidelity Master Mix

0.2ul of primer Xba-I-AarC-F

0.2ul of primer Repair-Pst-I-R

0.2ul of template DNA at 0.1ng

up to 10ul

(2) 5ul of I-5™ 2X High-Fidelity Master Mix

0.4ul of primer Repair-Pst-I-F

0.4ul of primer T7T-Pst-I-R

0.2ul of template DNA at 0.1ng

up to 10ul

PCR program

(1) 98C/2min

(2) 98C/10s

(3) 62C/10s

(4) 72C/20s

(5) Repeat cycle (steps 2 to 4, 34 times)

(6) 72C/2min

2. Splicing two fragment

Primers

gctctagatagctactagagtacacaggaagtagtagtgaccgaacgtatc Xba-I-AarC-F
tgactgcagcggcccgtactagtaaaaaaccctcaagaccgtagagg T7T-Pst-I-R

Diluted to 10 pmol/ul

Reaction system

(1) 25ul of I-5™ 2X High-Fidelity Master Mix

2ul of primer Xba-I-AarC-F

2ul of primer T7T-Pst-I-R

1ul of template DNA at 0.1ng (two fragments are equimolar)

up to 50ul

PCR program

- (1) 98C/2min
- (2) 98C/10s
- (3) 55C/10s
- (4) 72C/20s
- (5) Repeat cycle (steps 2 to 4, 34 times)
- (6) 72C/2min

HPLC

Sample

The bacterial cultures were centrifuged at 3000g for 10min.

Take the supernatant and centrifuged at 25000g for 10min.

Take the supernatant and filtered by a 0.22um filter.

Mobile phase

Methanol- 0.01mol/l K_2HPO_4 (3:97)

Phosphoric acid is used to regulate pH to 2.0

Chromatographic condition

C18 column(4.6×250mm, 5um)

Flow rate is 1ml/min

Injection volume is 60 ul

Detection wavelength is 210nm

Column temperature is 25C

Standard

Accurately weighed 1000mg acetic acid and use mobile phase adjust volume to 10ml

Dilute to 0.01mg/ml, 0.1mg/ml, 1mg/ml, 10mg/ml, 100mg/ml.

Use different concentration to draw standard curve.