


## DBCO-PEG4 - TAMRA

<b>Product No.:</b>	A131
<b>Product Name:</b>	DBCO-PEG <sub>4</sub> -TAMRA
<b>Alternative name:</b>	ADIBO-PEG <sub>4</sub> -TAMRA, DBCO-PEG <sub>4</sub> -Tetramethylrhodamine

**Chemical Structure:**



<b>Chemical Composition:</b>	C <sub>54</sub> H <sub>58</sub> N <sub>5</sub> O <sub>10</sub>
<b>Spectral Properties:</b>	Abs/Em = 545/565 nm
<b>Extinction Coefficient:</b>	92000 M <sup>-1</sup> cm <sup>-1</sup>
<b>Molecular Weight:</b>	937.07
<b>Appearance:</b>	Red solid
<b>Storage:</b>	Upon receipt store at -20°C. Product shipped at ambient temperature

DBCO-PEG<sub>4</sub>-TARMA is an azide-reactive fluorescent dye that is well suited for detection and labeling of chemically, enzymatically, or metabolically azide-modified biopolymers or peptides. The DBCO group reacts with an azide to produce a stable triazole (Figure 1), which is also referred to the Cu(I)-free or strain-promoted click reaction.




Figure 1. Reaction scheme of DBCO and azide

## Live Cell Labeling

1. Grow mammalian cells in an appropriate medium with an azide-derivatized metabolite (e.g., ManNAz) at 37°C in 5% CO<sub>2</sub>.
2. Wash the cells two times with D-PBS containing 1% FBS.
3. Prepare a 5 mM stock solution of DBCO-PEG<sub>4</sub>-TARMA in a water-miscible solvent such as DMSO or DMF by adding 0.427 mL of solvent to 2 mg vial or 1.07 mL to a 5 mg vial, and vortex to dissolve all solid.
4. Label the azide-modified cells at room temperature in the dark for 30-60 min with 5 to 30 μM of DBCO-PEG<sub>4</sub>-TARMA in D-PBS containing 1% FBS.
5. Wash the cells four times with D-PBS containing 1% FBS.
6. Fix the cells with 4% formaldehyde in D-PBS for 20 minutes at room temperature.
7. Wash the cells with D-PBS.
8. *Optional step:* Counterstain the cells for 15 minutes at room temperature with Hoechst 33342 in D-PBS.
9. Wash the cells two times with D-PBS.

10. Image the cells.

## Lysing Cells

**Do not use DTT, TCEP, or  $\beta$ -mercaptoethanol because they will reduce the azide.**

1. Prepare lysis buffer (100 mM Tris buffer (pH 8.0) containing 1% (w/v) SDS).
2. *Optional step:* Add protease and phosphatase inhibitors the lysis buffer.
3. Suspend cells in the lysis buffer (50  $\mu$ L lysis buffer per  $10^6$  cells) and heat to 75  $^{\circ}$ C. If using a 6-well plate you need 500  $\mu$ L lysis buffer per 100 mm dish and 200  $\mu$ L lysis buffer per well.
4. Sonicate the lysate briefly to shear DNA and reduce the viscosity of the solution.
5. Vortex the lysate for 5 minutes.
6. Centrifuge the cell lysate at 16,000 *g* at 4 $^{\circ}$ C for 10 minutes.
7. Transfer the supernatant to a clean tube and determine the protein concentration if required. Ideally, the protein concentration should be 1–2 mg/mL.
8. Prepare 1 M solution of iodoacetamide by adding 3 mL of DMSO to IAA labeled vial (provided with a lysis labeling kit).
9. Block cysteine thiols in lysate by addition of iodoacetamide stock solution to a final concentration of 15 mM, agitate mildly for 30 min.
10. Prepare a 5 mM stock solution of DBCO-PEG<sub>4</sub>-TARMA by adding 0.427 mL of DBCO to 2 mg vial or 1.06 mL to 5 mg vial.
11. Label the cell lysate by addition of DBCO-PEG<sub>4</sub>-TARMA to a final concentration of 20  $\mu$ M. Protect from light and agitate mildly for 30 min at room temperature.
12. Prepare a 50 mM stock solution of stop buffer by adding 3 mL of water to Stop Reagent labeled vial (provided with a lysis labeling kit).
13. Stop reaction by addition of stop buffer to a final concentration of 100  $\mu$ M, agitate briefly for 20 min.
14. Load  $\sim$ 10  $\mu$ g of protein on 12% Tris-Tricine SDS-PAGE gel.
15. Image the gel by fluorescence scanning with detection for TAMRA (TMR), or AF 546.
16. Stain the gel with Coomassie-stain according to the manufacture's protocol, and image cells.

## References

1. (a) Baskin, J. M., *et al.* (2007). Copper-free click chemistry for dynamic *in vivo* imaging. *PNAS.*, **104**(43):16793-7. (b) Ning, X., *et al.* Visualizing metabolically labeled glycoconjugates of living cells by copper-free and fast Huisgen cycloaddition. *Anweg. Chem. Int. Ed.*, **47**:2253-5.

2. (a) Rubino F. A., *et al.* (2012). Chemoselective Modification of Viral Surfaces via Bioorthogonal Click Chemistry. *J. Vis. Exp.*, **66**:4246. (b) Yao J. Z., *et al.* (2012). Fluorophore Targeting to Cellular Proteins via Enzyme-Mediated Azide Ligation and Strain-Promoted Cycloaddition. *J. Am. Chem. Soc.*, **134**:3720–3728.
3. (a) Dieterich, D.C., *et al.* (2007). Labeling, detection and identification of newly synthesized proteomes with bioorthogonal non-canonical amino acid tagging. *Nature Protocols*, **2**(3): 532-40. (b) Best, M.D., *et al.* (2009). Click chemistry and bioorthogonal reactions: unprecedented selectivity in the labeling of biological molecules. *Biochemistry*, **48**(28):6571-84. (c) Ngo J. T., *et al.* (2012). State-selective metabolic labeling of cellular proteins. *ACS Chemical Biology*, **7**(8):1326-30.