

GROWTH CONDITION, STORAGE, MEDIA & ANTIBIOTICS FOR *BACILLUS SUBTILIS*

Growth conditions:

B. subtilis growth best at 37°C and has a doubling time of 30 min. Always use flasks that comprise at least 5x times the volume of media used, and always use lids that are able to allow air passage. *B. subtilis* grows strictly aerobic!

Storage:

We use three types of storage conditions for *B. subtilis* cells.

1. **Glycerol-stocks** for long-term storage and back-up: Freeze two independent clones of your *B. subtilis* strain at -80°C with glycerol (1.3 ml fresh overnight culture + 50% glycerol).
2. **DSM-Plate:** Long-term storage for frequent used strains. Strike out your cells on a DSM plate. *B. subtilis* is able to produce spores, that are very resistant to all kind of environmental conditions. So storage of spores is a reliable way of storage.
3. **LB-plates:** For everyday use you can leave your *B. subtilis* strains on LB-plates at room temperature. **Note:** *B. subtilis* is cold sensitive, never store *B. subtilis* cells at 4°C.

Media

Luria-Bertani (LB) broth:

Tryptone	10 g
Yeast extract	5 g
NaCl	10 g
H ₂ O (dest)	ad 1.000 ml

- for LB plates: add 15 g/l of agar
 - important: cool down the agar solution to 50°C before adding antibiotics

Difco Sporulation Medium (DSM):

Nutrient Broth	8 g
----------------	-----

KCl	1 g
MgSO ₄ (1 M)	1 ml
MnCl ₂ (10 mM)	1 ml
H ₂ O (bidest)	ad 1.000 ml
Add after autoclave:	
CaCl ₂ (1 M)	0,5 ml
FeSO ₄ (1 mM)	1 ml

- for DSM plates: add 15 g/l of agar
 - important: cool down the agar solution to 50°C before adding antibiotics

Starch plates:

Nutrient Broth (Difco)	7,5 g
Starch	5 g
Agar	15 g
H ₂ O (dest)	ad 1.000 ml

Chemical defined medium (CSE):

1×C-Salts	100 ml
Tryptophan (5 mg/ml)	1 ml
Ammoniumeisencitrat (2,2 mg/ml)	1 ml
III'-Salts	1 ml
Potassium glutamate (40%)	2 ml
Sodium succinate (30%)	2 ml

5×C-Salts (1 l)

KH ₂ PO ₄	20 g
K ₂ HPO ₄ × 3 H ₂ O	80 g

$(\text{NH}_4)_2 \text{SO}_4$	16,5 g
-------------------------------	--------

III'-Salts (1 l)

$\text{MnSO}_4 \times 4 \text{H}_2\text{O}$	0,232 g
---	---------

$\text{MgSO}_4 \times 7 \text{H}_2\text{O}$	12,3 g
---	--------

- autoclave (or filtrate) each component separately and put them together freshly before starting your experiment
- Optionally: addition of media additives, for example pyruvate (0.5% final concentration) or glucose (1% final concentration)

Antibiotics

- Indicated are 1.000-times stock solutions
- Dissolve in the specific solvent and filtrate by using 0.2 µm filters
- Store at -20°C

Strain	Antibiotic	Concentration	Dissolve in	Color code
<i>B. subtilis</i>	Kanamycin	10 mg/ml	H ₂ O	Black (one bar)
	Chloramphenicol	5 mg/ml	70% ethanol	Blue
	MLS selection:			Red
	Erythromycin	1 mg/ml	70% ethanol	
	Lincomycin	25 mg/ml	H ₂ O	
	Spectinomycin	100 mg/ml	H ₂ O	Purple
	Tetracyclin	10 mg/ml	H ₂ O	Orange
	Bacitracin	50 mg/ml	H ₂ O	-
<i>E. coli</i>	Ampicillin	100 mg/ml	H ₂ O	Green
	Kanamycin	50 mg/ml	H ₂ O	Black (two bars)
	Streptomycin	20 mg/ml		
<i>Rhodobacter</i>	Tetracyclin	1 mg/ml	H ₂ O	Orange
	Spectinomycin	25 mg/ml	H ₂ O	Purple
	Kanamycin	25 mg/ml	H ₂ O	Black

Further substances

- IPTG: 1 M (dissolve in H₂O; store at -20°C)
- X-Gal: 100 mg/ml (dissolve in DMF; color code: blue X; store at -20°C)
- BCIP (5-Brom-4-chlor-3-indoxylphosphat-p-toluidinsalz): 50 mg/ml (dissolve in DMF; store at -20°C)