

# Preparation of Protoplasts

1. For *Cochliobolus*, inoculate 100 mL complete medium (CM) in a 500 mL flask with  $10^6$ - $10^7$  conidia. Alternatively, mycelium from solid or liquid medium may be blended for 5-10 seconds in a sterile blender cup and then used as inoculums. For other fungi, grow as usual on medium known to produce conidia. 150  $\mu$ g/mL ampicillin may be used to prevent bacterial contamination.
2. Grow 12-18 hours at room temperature (25-28°C) with shaking at 150-250 rpm.
3. Transfer the culture to 40 mL centrifuge tubes and centrifuge at 5000 rpm for 5 minutes at 4°C in a SS34 rotor (Sorvall) or its equivalent. Discard the supernatant and try not to disturb the pellet.
4. Scoop 0.5-0.8 g portions of the pellet and resuspend well in 10 mL of enzyme-osmoticum in a 50 mL Erlenmeyer flask (use 80 mL total of enzyme-osmoticum, 8 x 10 mL).
5. Place flasks at 30-32°C and shake gently (50 rpm) for 1.5-2 hr.  $10^8$ - $10^9$  protoplasts should be released in this time (count with a hemacytometer).
6. Separate the protoplasts from the intact mycelium and cell wall debris by filtering the protoplast suspension through 4 layers of sterile cheesecloth, then through sterile nylon fabric with a 20-25  $\mu$ m pore size.
7. Pellet the protoplasts by centrifuging the filtrate from step 6 in 40 mL screw-capped centrifuge tubes at 5000 rpm for 5-6 min in a SS34 rotor at 4°C. Discard the supernatant very carefully trying not to disturb the protoplast pellet. The enzyme-osmoticum can be re-used twice if stored at -20°C.
8. Resuspend the pellets gently in a total of 10 mL 0.7 M NaCl; combine into one tube. Recentrifuge as in step 7. Discard the supernatant.
9. Wash the pellet 3 times with 10 mL STC, pelleting the cells between washes by centrifuging as in step 7. Resuspend protoplasts in 500  $\mu$ L STC after final wash.
10. Count protoplasts and adjust their concentration to approximately  $10^8$ /mL with STC.