

# Preparation of *Bacillus subtilis* competent cell


---

## Materials

- Sterile 1.5 mL microtubes
- LB medium
- T base – Autoclave
  - $(\text{NH}_4)\text{SO}_4$  2g
  - $\text{K}_2\text{HPO}_4 \cdot 3\text{H}_2\text{O}$  18.3g
  - $\text{KH}_2\text{PO}_4$  6g
  - Trisodium citrate .  $2\text{H}_2\text{O}$  1g
- SpC – Made fresh on day
  - T base 20ml
  - 50%(w/w) glucose 0.2ml
  - 1.2%(w/w)  $\text{MgSO}_4 \cdot 3\text{H}_2\text{O}$  0.3ml
  - 10%(w/w) Bacto yeast extract 0.4ml
  - 1%(w/w) tryptone 0.5ml
- SpII – Made fresh on day
  - T base 200ml
  - 50%(w /w) glucose 2ml
  - 1.2%(w/w)  $\text{MgSO}_4 \cdot 3\text{H}_2\text{O}$  14ml
  - 10% Bacto yeast extract 2ml
  - 1%(w/v) tryptone 2ml
  - 0.1 M  $\text{CaCl}_2$  1ml

## Apparatus

- Centrifuge
- Spectrofotometer
- Shaker
- Flow chamber
- Autoclave



## Method

1. Streak out the strain to be made competent on an LB agar plate as a large patch and incubate overnight at 30°C;
2. The following morning scrape the cell growth off the plate and use to inoculate fresh, pre-warmed, SpC medium (20 mL) to give an OD<sub>600</sub> reading of about 0.5;
3. Incubate the culture at 37 °C with vigorous aeration and measure OD periodically (OD<sub>600</sub>) to assess cell growth;
4. When the rate of cell growth is seen to depart from exponential (i.e. no significant change in cell density over 20-30 min) inoculate 200 mL of pre-warmed, SpII medium with 2 mL of stationary-phase culture and continue incubation at 37 °C with slower aeration;
5. After 90 min incubation, pellet the cells by centrifugation (4,000 g, 25min) at room temperature;
6. Carefully decant the supernatant into a sterile container and save;
7. Resuspend the cell pellet in 2,25 mL of the saved supernatant, by mixing gently.