## Cultivation protocol of Clostridium stercorarium

Version: v1.1 Release: 3rd Sep.

We received the freeze dried *Clostridium stercorariumculture* (DSM No.8532) delivered from DSMZ in ampoules.

To cultivate *C. stercorarium*, culture medium should be prepared first. The recipe for *C. stercorarium* is shown as below:

$KH_2PO_4$	<b>0.50</b> g
K <sub>2</sub> HPO <sub>4</sub> ·3 H <sub>2</sub> O	1.00 g
Urea	2.00 g
MgCl <sub>2</sub> ·6 H <sub>2</sub> O	<b>0.50</b> g
CaCl <sub>2</sub> ·2 H <sub>2</sub> O	<b>0.05</b> g
FeSO <sub>4</sub> ·7 H <sub>2</sub> O	1.25 mg
Morpholinopropane sulfonic acid	10.00 g
Resazurin	1.00 mg
Yeast extract	6.00 g
Cellobiose	5.00 g
Cysteine-HCl·H <sub>2</sub> O	1.00 g
Distilled water	1000.00 ml
Adjust pH	I to 7.2

## Decolor Four iGEM Tongji 2014

The ingredients except cellobiose were mixed well by swirling and autoclaved at 121°C for 25 min while cellobiose was sterilized separately by filtration. The liquid medium (or broth) could be used for growing bacteria in test tubes.

Broth with about 1.5% agar added was poured into Petri dishes and let each plate cool until totally solid. The plates could be used to isolate pure culture in the following step.(Figure 1)



Figure 1




Figure 2

When the medium was ready, we opened the ampoules carefully and rehydrated the dried cultures near a flame to keep it under sterile conditions. (Figure 2)

Then, **0.5ml** culture medium was added in the vial to allow the pellet to rehydrate for **30 minutes**. The whole content was mixed gently and transferred to a **15ml tube with 5ml liquid medium**. The bacteria were cultivated at **60°C overnight** in an anaerobic bag.

The following day, a small sample from the *Clostridium* stercorarium inoculums was spread on the surface of the plate to separate single colonies by using spread plate method (no dilution). The plates were put into an anaerobic bag and cultured at 60°C for 2 days.

However, we did not get single colony because of the high density of colonies at first. Then, 10-3, 10-4 and 10-5 dilutions of the original broth were respectively spread on the plates. Finally, single colonies were successfully isolated from these plates and could be applied for further study.