

Cultivation and transformation of *Synechocystis* strains.

Strain	Genotype	Reference
WT <i>Synechocystis</i> sp. PCC 6803		

Growth media: BG-11
WT is glucose tolerant

Preparation of liquid BG-11 media [1]:

Add 1.5 g NaNO₃ to 1 L of ddH₂O.

Optional: Add 5.95 g/L HEPES

Add 1 ml of each filter sterilized stock solution below.

Set pH to 7.5-7.8.

Autoclave.

Stock 1 - 50 ml ddH ₂ O	
K ₂ HPO ₄	1.53 g
Stock 2 - 50 ml ddH ₂ O	
MgSO ₄ *7 H ₂ O	3.75 g
Stock 3 - 50 ml ddH ₂ O	
CaCl ₂ *2 H ₂ O	1.8 g
Citric acid	0.3 g
Stock 4 - 50 ml ddH ₂ O	
Ferric ammonium citrate	0.3 g
EDTA (disodium magnesium salt)	0.05 g
Stock 5 - 50 ml ddH ₂ O	
Na ₂ CO ₃	1 g
Stock 6 - 1000 ml ddH ₂ O	
H ₃ BO ₃	2.86
MnCl ₂ *4 H ₂ O	1.81
ZnSO ₄ *7 H ₂ O	0.222
Na ₂ MoO ₄ *2 H ₂ O	0.39
CuSO ₄ *5 H ₂ O	0.079
Co(NO ₃) ₂ *6 H ₂ O	0.0494

Optional: Prepare a stock solution of 0.5 M NaHCO₃ in BG-11, e.g. 100 ml:

4.2 g NaHCO₃

0.595 g HEPES

100 ml ddH₂O

100 ul of each of the six BG-11 stock solutions

Adjust pH to around 7.8 with NaOH*

Do NOT autoclave, filter through 0.2 μm filter

*pH typically ends up at about 7.7-7.8 when adding this much NaHCO_3 , so further pH adjustment is usually not needed.

Then you just mix 1/10 volume of stock with 9/10 volume of your normally prepared BG-11 supplemented with HEPES. You can always double check pH with pH paper in the media you have used so far (maybe for NaHCO_3 media in particular) if you think pH might have been a problem.

Optional: A buffer should be added if cultures are grown in higher CO_2 than atmospheric levels, or if grown with up to 50 mM NaHCO_3 supplemented. Use e.g. 25 mM HEPES (5.95 g/L) for this purpose. *Synechocystis* can also be grown on glucose (e.g. 20 mM), which is necessary at dark, anoxic conditions.

Preparation of BG-11 agar plates:

For 200 ml

Prepare two flasks to be autoclaved separately:

1

- 3 g agar (standard bacteriological agar)
- 0.3 g NaNO_3
- 100 ml ddH₂O

2 – set pH to 7.5-7.8

- 0.6 g sodium thiosulfate
- 100 ml ddH₂O
- 200 μl of each stock solution
- optional: 1.19 g HEPES

Mix the two after autoclaving and add antibiotics as appropriate at 50 °C.

Alternatively, the two flasks can be stored separately in room temp for weeks/months. When agar plates are needed, just microwave the flask with agar (the flask should preferably not be filled to more than half) until all agar has melted. Mix 1:1 with room temp liquid phase (flask 2). The mix will then be at approximately 50 °C and antibiotics can be added directly.

Use approximately 20 ml/plate (larger volumes, e.g. 30 ml, can be used for long-term storage).

Plates *can* be sealed with one layer of parafilm to keep a moister environment. This also prevents the plate from drying out too fast.

Grow plates and liquid cultures at 20-30 °C (faster growth at 28-30 °C). Light intensities of about 50 $\mu\text{mol photons m}^{-2} \text{s}^{-1}$ are commonly used, but cultures can grow well both at higher and lower intensities. Very dilute cultures or recently plated transformants may be a bit sensitive to *strong* light and can be protected by a thin layer of paper if necessary.

Doubling time depends on the growth conditions, but is typically 10-12 h at around 50 $\mu\text{mol photons m}^{-2} \text{s}^{-1}$ and 1% CO_2 or in media supplemented with glucose or NaHCO_3 .

Natural transformation

Inoculate a culture to OD₇₃₀ 0.1 or less. Grow until approximately OD₇₃₀ 0.7. The plasmid used should be non-replicative in *Synechocystis* and have two homologous regions (ca 1000 bp each) for integration into the genome.

Protocol:

- Centrifuge 10 mL of healthy cells at an O.D_{730nm} of ~ 0.7 at 4000xg for 6 min
- Wash pellet twice in BG-11
- Resuspend the cell pellet in 200 µL of BG-11 in a 2 ml eppendorf tube
- Add 1 µg of target plasmid to the cells and incubate the tube in light for 5 hours at 30 °C.
- Add 4 ml BG-11 and grow in a flask 18-24 h
- Plate cells on BG-11 plates with appropriate antibiotics

References

- [1]: Rippka R, Deruelles J, Waterbury JB, Herdman M, Stanier RY. Generic assignments, strain histories and properties of pure cultures of cyanobacteria. Microbiology. 1979 Mar 1;111(1):1-61.