

## iGEM 2017 – Microbiology – BMB – SDU

**Title:** LA plates with antibiotic

**SOP number:** SOP01

**Version number:** 03

**Date issued:** 2012.11.26

**Review date:** 2017.09.10

**Original by:** TKR

**Modified by:** FN, EG

### 1. **Purpose**

To make LA plates with selection

### 2. **Area of application**

This procedure is valid for all antibiotics

### 3. Apparatus and equipment

Apparatus/equipment	Location (Room number)	Check points	Criteria for approval/rejection
Refrigerator	Laboratory (class 1)	•	121°C
Small autoclave	Laboratory (class 1) – V18-405-0	•	
Pipette (p1000, 100)		•	
Measuring pitcher 1 L	Glass hallway	•	

### 4. Materials and reagents – their shelf life and risk labelling

Name	Components (Concentrations)	Manufacturer / Cat. #	Room	Safety considerations
LA	1% Tryptone 1% NaCl 0.5% Yeast extract 1.5% agar	Oxoid Sigma-Aldrich Merck Difco agar from BD	Anne-Mette Or V18-405-0	
Appropriate antibiotic	Manufacture/cat#, stock solution conc.; Stock solution date of creation, dilution in medium		V18-405-0	
Plastic Petri dishes		Contact lab-manager (Sarstedt)	BMB storage	
Green Pipette tips		Contact lab-manager	Micro storage	
Blue Pipette tips		Contact lab-manager	Micro storage	
Plastic bag	LDPE standardpose	CC&CO Living Office Lev varenr: 281615		

## 5. QC – Quality Control

Calculation of antibiotic concentration

*Antibiotic stock concentration*  $\times$  *Required volumen of antibiotica*

$=$  *Desired concentration of antibiotica*  $\times$  *Total volumen of media*

$\Rightarrow$

*Required volume of antibiotica*  $= \frac{\text{Desired concentration of antibiotica} \times \text{Total volumen of media}}{\text{Antibiotic stock concentration}}$

The temperature of the LA media may not be more than 60°C when the antibiotic is added

It is possible to leave the autoclaved LA media at 58°C up till 24 hours before adding antibiotic and pouring into plates.

Note that some antibiotic are decomposed over time, specially ampicillin

## 6. List of other SOPs relevant to this SOP

iGEM2017\_SOP14\_v02\_EG\_Table\_Autoclave

## 7. Environmental conditions required

## 8. Procedure

Mix components of LA in glass bottle and autoclave (amount is labeled on)

- 8.1 Add the appropriated amount of antibiotic when the temperature is between 50-60 °C, see appendix
- 8.2 Fill up the required amount of petri dishes (between 15-25 ml LA in each)
- 8.3 Dry on table for at least 2 hours, can stand ON
- 8.4 Note type and concentration of antibiotic and the date at the bottom of the dishes
- 8.5 Place in bag in the refrigerator until use

## 9. Waste handling

Chemical name	Concentration	Type of waste (C, Z...)	Remarks
LA media rest		GMO waste	

## 10. Time consumption

- Total-time 3½ hours if the dishes are dried 3 hours
- Hands-on-time ½ an hour

## 11. Scheme of development

Date / Initials	Version No.	Description of changes
12.11.26 / TKR	01	The SOP has been written
13.05.22 / MM & TK	01	The SOP has been approved
17.04.28 / FN	02	The SOP has been modified
17.05.01 / EG	02	The SOP has been approved
17.09.10 / EG	03	The SOP has been modified

## 12. Appendices

### - Addition of antibiotic

- Chloramphenicol
  - Add 1  $\mu$ L of 30 mg/mL chloramphenicol per mL LA/LB to achieve 30  $\mu$ g/mL.
- Ampicillin
  - Add 1.0  $\mu$ L of 100 mg/mL ampicillin per mL LA/LB to achieve 100  $\mu$ g/mL.
- Kanamycin
  - Add 1.0  $\mu$ L of 50 mg/mL kanamycin per mL LA/LB to achieve 50  $\mu$ g/mL.