

## PureYield™ Plasmid Miniprep System from Promega

- ◆ Prepare Lysate
  - ◇ Add 600  $\mu\text{L}$  of bacterial culture to a 1.5 mL microcentrifuge tube.
  - ◇ Add 100  $\mu\text{L}$  of **Cell Lysis Buffer (Blue)**, and mix by inverting the tube 6 times.
  - ◇ Add 350  $\mu\text{L}$  of cold **Neutralization Solution**, and mix thoroughly by inverting.
  - ◇ Centrifuge at maximum speed in a microcentrifuge for 3 minutes.
  - ◇ Transfer the supernatant to a **PureYield™ Minicolumn** without disturbing the cell debris pellet.
  - ◇ Place the minicolumn into a Collection Tube, and centrifuge at maximum speed in a microcentrifuge for 15 seconds.
  - ◇ Discard the flowthrough, and place the minicolumn into the same Collection Tube.
- ◆ Wash
  - ◇ Add 200  $\mu\text{L}$  of **Endotoxin Removal Wash (ERB)** to the minicolumn. Centrifuge at maximum speed in a microcentrifuge for 15 seconds.
  - ◇ Add 400  $\mu\text{L}$  of **Column Wash Solution (CWC)** to the minicolumn. Centrifuge at maximum speed in a microcentrifuge for 30 seconds.
- ◆ Elute
  - ◇ Transfer the minicolumn to a clean 1.5 mL microcentrifuge tube, then add 30  $\mu\text{L}$  of **Elution Buffer** or nuclease-free water directly to the minicolumn matrix. Let stand for 1 minute at room temperature.
  - ◇ Centrifuge for 15 seconds to elute the plasmid DNA. Cap the microcentrifuge tube, and store eluted plasmid DNA at  $-20\text{ }^{\circ}\text{C}$ .

From: [Promega](#)