

IDT DNA.com Directions for miniGENE

Plasmid Resuspension

1. Centrifuge tube before opening
2. Resuspend DNA in 80.0 μL of IDTE buffer
 - should reach a concentration of 50.0 ng/ μL
 (Total of 4000 ng in tube (4 μg))
3. Incubate @ room temp for 30 minutes.
4. Centrifuge for 2 minutes.

o Dilution for transformation:

- o add 2.0 μL DNA to 999.0 μL diH₂O
 \rightarrow will provide 0.1 ng/ μL concentration.
- o use 1-2 μL of dilution for each transformation.

o Dilution for PCR concentration: 10 ng/ μL

use 1.0 μL of 50 ng/ μL DNA, mix w/ 4.0 μL diH₂O.
 = 5.0 μL 10 ng/ μL

IDT Primer Dilutions:

"Mamba EcoRI Fwd" 3.0 nmol

Resuspend w/ 300.0 μL \rightarrow 10.0 μM concentration.

"Mamba NotI Rev" 2.9 nmol

Resuspend w/ 290.0 μL \rightarrow 10.0 μM concentration

do aliquots 20 μL of 10 μM conc

2 μL stock + 18 μL diH₂O ✓

both ✓

Neb Q5 High Fidelity 2x master mix PCR

Purpose: Amplify mamba cds sequence from JDT minichrom vector while also adding EcoRI and NotI restriction sites.

Reaction mixture Total 50.0 μ L

25.0 μ L Q5 High Fidelity 2x master mix
 2.5 μ L mamba EcoRI Fwd 10 μ M
 2.5 μ L mamba NotI Rev 10 μ M
 1.0 μ L DNA template (10 ng/ μ L)
 + 19.0 μ L Nuclease free dH₂O
 50.0 μ L total volume.

⊖ Negative control rxn

25.0 μ L Q5 High Fidelity 2x master mix
 2.5 μ L mamba EcoRI Fwd
 2.5 μ L mamba NotI Rev
 + 20.0 μ L Nuclease Free dH₂O
 50.0 μ L Total volume.

Amplification Protocol

• Thermocycler Settings.

98°C (30) seconds

25-35 cycles { 98°C 5-10 seconds
 58°C 10-30 seconds
 72°C 20-30 seconds } ~ 43 mins

72°C 5.0
 4°C 2.0 minutes } 1 hr 21 mins