

## Growth inhibition assay

### Aim

- To establish whether the transformation products (TPs) generated from the reactions with laccase, SMX with or without the mediator ABTS are more toxic than pure SMX
- To test if the bacteria *Bacillus subtilis* (*B. subtilis*) and *Escherichia coli* (*E. coli*) can be used

### Materials

- *B. subtilis* strain ISW1214, *E. coli* strain TOP10
- LB media (see separate protocol for making LB media)
- TPs from reactions (see separate protocol for making reactions)
- 96 well plate
- FLUOstar Omega analyser (BMG LABTECH)

### Procedure (from plates)

- Streak *B. subtilis* and *E. coli* on separate LB plates (without antibiotics) and incubate the plate overnight
- Pick cells and inoculate in 20 mL of LB media (two separate shaker flasks)
- Let the cells grow overnight at 37 °C, between 150-200 rpm
- Dilute the *B. subtilis* and the *E. coli* to an OD600 of 0.1/ mL LB media
- Pipette 185 µL of cells into wells, and 15 µL of transformation products or pure SMX on the cells
- Pipette them into a 96 well plate and measure OD600 in a plate reader (FLUOstar Omega, BMG LABTECH) every 5 minutes for 7 hours

### Notes

Make all samples in triplicates (n=3). Use a LB medium triplicate as a blank. Use pure SMX and cells as control. Pipette pure cells into a triplicate to compare with “ultimate” growth.

After culturing the cells overnight on the plate, it is possible to store the plates in 4 °C for further use.

### References

Becker, D., Varela Della Giustina, S., Rodriguez-Mozaz, S., Schoevaart, R., Barceló, D., de Cazes, M., ...

## Growth inhibition assay

Wagner, M. (2016). Removal of antibiotics in wastewater by enzymatic treatment with fungal laccase – Degradation of compounds does not always eliminate toxicity. *Bioresource Technology*, 219, 500–509. doi.org/10.1016/j.biortech.2016.08.004

iGEM Bonn, iGEM F. (2016). *What is Bacillus subtilis?* Retrieved from 2016.igem.org/wiki/images/7/74/T--UBonn\_HBRS--How-To-Bacillus-Subtilis.pdf

iGEM Toulouse. (2016). *Study of the optimum growth conditions of Bacillus subtilis (strain WT 168), and Pseudomonas fluorescens (strain SBW25)*. Retrieved from 2016.igem.org/wiki/images/e/e0/T--Imperial\_College--ToulouseGrowthStudy.pdf