

Red Blood Cell Preparation and Hemagglutination Assay

Purpose: This protocol describes the preparation of RBCs for storage and use in Hemagglutination assays, which are used to determine the La Sota B1 lentogenic Newcastle Disease virion hemagglutinin-neuraminidase titer relative to virion stock positive controls. The protocol was adapted from McGinnes et al. (2006) and “Detection of Hemagglutinating Viruses” (n.d.).

Materials:

- Whole Blood Cells
- 50mL conical tube
- Multichannel micropipette, micropipette, and tips
- PBS with 2mM Penicillin/Streptomycin
- SV3
- Alsever's Solution
- 96 round bottom well plate
- Microscope slide
- Virion stock
- Centrifuge

Procedure:

1. RBC preparation
 - a. Obtain 25mL whole blood and add 25mL cold Alsever's Solution for anti-coagulation in a 50mL conical tube and keep on ice
 - b. Centrifuge whole blood at 500 RCF for 10min
 - c. Aspirate blood plasma, buffy layer, and top erythrocytes
 - d. Wash 3 times by resuspending in PBS with 2mM Penicillin/Streptomycin to double the pellet volume, centrifuging at 500 RCF, and aspirating the supernatant
 - e. Resuspend the pellet in PBS with 2mM Penicillin/Streptomycin or SV3 for a final concentration of 10% pellet volume per total volume, and store at 2-7°C for 1 week or 42 days
2. HA Titer
 - a. Thaw your virion stock and samples on ice
 - b. Use a multichannel micropipette to add 50µL of cold PBS to each row on a 96 round bottom well plate for each sample in duplicates
 - c. Add 50µL of each sample to column 1 of their respective duplicate rows
 - i. Don't add anything to the negative control rows and add virion stocks to the positive control rows
 - d. Use a multichannel micropipette to make a 2-fold serial dilution on every row except for the negative control
 - e. Add 50µL of 1% RBC solution to each well and agitate the plate to mix
 - f. Cover plate in sterile rap and incubate for 1hr at 2-7°C

- g. Observe HA titer of each sample from the inverse of the last serial dilution that indicates RBC agglutination
 - i. Agglutinated RBCs will appear cloudy or be followed by cloudy wells in further dilutions as a result of the hook effect
 - ii. Non-agglutinated RBCs will have a red button at the bottom of the well, where RBCs settled
- 3. HA Test
 - a. Thaw your virion stock and samples on ice
 - b. Add a drop of blood to a microscope slide for the positive control (Virion Stock), the negative control (PBS), and each sample
 - c. Add a drop of undiluted samples and controls to their respective drops of blood
 - d. Agitate the slide to mix for 1min
 - e. Observe agglutination under the microscope
 - i. Positive: chunky lattice structure forms
 - ii. Negative: blood is closely packed and homogeneous

References

McGinnes, L.W., Pantua, H., Reitter, J., and Morrison, T.G. (2006). Newcastle Disease

Virus: Propagation, Quantification, and Storage. In Current Protocols in Microbiology, R. Coico, A. McBride, J.M. Quarles, B. Stevenson, and R.K. Taylor, eds. (Hoboken, NJ, USA: John Wiley & Sons, Inc.), pp. 15F.2.1-15F.2.18.

"Detection of Hemagglutinating Viruses." Edited by Linda S. Snively , *United States Department of Agriculture Center for Veterinary Biologics Testing Protocol*, United States Department of Agriculture Animal and Plant Health Inspection Service, 9 Jan. 2018, www.aphis.usda.gov/animal_health/vet_biologics/publications/VIRPRO0096.pdf.