

Cleaning Glass Tubes

Overview

This protocol is for the clean up for single use glass tubes and their plastic caps.

Materials

- 10% bleach can be made using the stock solution of bleach in the cleaning supplies cabinet
- Dirty glass tubes

Procedure

1. Using the bleach solution bottles, slightly fill the glass tubes and let sit 1-3 minutes
2. Dump out the glass tubes down the sink and dispose of the glass tubes in the glass waste container
3. Place the plastic tops in a waste container and submerge them in 10% bleach solution
4. Let sit for 5 minutes
5. Rinse thoroughly with water
6. Rinse again with Millipore water
7. Let dry
8. Place back in plastic cap bin