

## **Refilling Bleach Bottles**

### **Overview**

This protocol covers refilling the lab bleach squirt bottles.

### **Materials**

- DI water
- Bleach Bottles
- Bleach (Sodium Hypochlorite)
  - o Located under the back room sink.

### **Procedure**

1. Bottles are to be filled with ~10% v/v bleach.
  - o Be careful. Don't ruin your clothes.
2. Measure out roughly 10% of the bottle's volume in bleach (roughly slightly below the first textured line).
3. Add DI water to fill line.
  - o Go slowly and be careful. The DI comes out very fast. Pouring water too fast will get your clothes bleached, and will lead to bubbles forming, which causes refilling to take longer.