




# 3 Parallel Droplet Generator and Mixer System

## Table of Contents


<b>Table of Contents</b>	<b>1</b>
<b>Overview</b>	<b>2</b>
<b>Chip Design</b>	<b>3</b>
<b>Milling Instructions</b>	<b>3</b>
<b>Testing Protocol</b>	<b>5</b>
Flow Layer Set Up	5
Testing the Chip	6
Cleaning the Chip	6

# Overview




This chip is designed to be used with a TERRA Adapter. It consists of three identical parallel systems, each containing a mixer and a droplet generator. It can run one, two, or all three of these systems at any one time.

# Chip Design



Flow Layer

## Milling Instructions




Flow Layer

## Notes:

- This chip should be milled on thick polycarbonate ( $5.00\text{mm} < Z_{\text{Polycarbonate}}$ )
- This chip should be used with thin PDMS ( $0.24\text{mm} < Z_{\text{PDMS}} < 0.26\text{mm}$ )

All the required SVGs for milling this chip are provided in the ZIP file. The layer, depth, and tool required for each SVG are listed in the file name. Below is the key describing how to read an SVG file name.




Mill the layers in the order they are listed with the correct depths and using the correct tool

Flow Layer	
Order	Layer Name
1	F_400_32
2	F_400_64
3	F_400_100
3	F_THRU_8
4	Border



# Testing Protocol

## Flow Layer Set Up



Inputs		
Name	Liquid	Flow Rate
A, B, C	Mineral Oil	4.0 mL/hour
D, E	Mixer 1 Inputs	2.0 mL/hour
F, G	Mixer 2 Inputs	2.0 mL/hour
H, I	Mixer 3 Inputs	2.0 mL/hour

<b>Outputs</b>	
<b>Name</b>	<b>Liquid</b>
<b>1</b>	Mixed Droplets from Mixer 1 in Mineral Oil
<b>2</b>	Mixed Droplets from Mixer 2 in Mineral Oil
<b>3</b>	Mixed Droplets from Mixer 3 in Mineral Oil

### Testing the Chip: Set Up and Protocol

1. Attach all the syringes to their correct inputs and set each to the desired flow rate indicated in the above table
2. Attach the output tubing to the chip's output
3. Start flowing all inputs through the chip
4. Let stabilize
5. Attach the outputs to the inputs of the TERRA Adapter (see its documentation for more information)

### Cleaning the Chip

1. Carefully disconnect tubing and dispose of all liquid waste
2. Disconnect all syringes
3. Clean the chip following the MARS protocols
4. Store the chip following the MARS protocols